
Instructions for Use

COLD RUSH®
Cold Therapy System

3

EN | Instructions for Use




5

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

16

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28

EN – Caution: Össur products and components are designed and tested according to the applicable official standards or an in-house defined standard when no official standard applies. Compatibility and compliance with these standard is achieved only when Össur products are used with other recommended Össur components. If an unusual movement or product wear is detected in a structural part of a device at any time, the patient should be instructed to immediately discontinue use of the device and consult his/her clinical specialist. This product has been designed and tested based on single patient usage. This device should NOT be used by multiple patients. If any problems occur with the use of this product, immediately contact your medical professional.

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Consult Instructions for use

Meets IEC type BF leakage current requirements

IP21

Protected against solid foreign objects of 12.5mmØ and greater. Protected against vertically falling water drops.

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ENGLISH

SAFETY PRECAUTIONS

WARNING: The Cold Rush® cold therapy system can be cold enough to seriously injure skin. Read the instructions which are located inside the unit, the instructions for the pad, and the labels on the device carefully before operating the Cold Rush cold therapy system.

WARNING: RX ONLY.

The device should only be used with a prescription from a physician that includes the following treatment information: (1) the duration of the entire cold therapy protocol, (2) the length of individual cold therapy sessions, as well as the length of breaks from cold therapy treatment, (3) setting of flow control knob, (4) instructions about how to inspect the skin, and (5) frequency of skin checks.

WARNING: SKIN INSPECTION AND ADVERSE REACTIONS.

The patient or his healthcare provider must check the skin at least every 1-2 hours during use of this device, regardless of whether the patient is asleep or awake. The patient must discontinue the use of the Cold Rush Cold Therapy Unit if he experiences any adverse reaction to the treatment area such as: increase in pain, burning, blistering, itching, increased swelling, discoloration of the skin, increased redness, welts or any other changes in the appearance of the skin.

WARNING: NO DIRECT CONTACT WITH SKIN:

INSULATION BARRIER REQUIRED (not included).


The Cold Rush Pad should never make direct contact with the skin. An insulation barrier between the pad and the skin should be used at all times without exception in accordance with instruction from the physician so as to prevent the skin from becoming too cold. The Cold Rush pad operating temperature is too cold to be placed directly on the skin.

- Unit is not suitable for use in the presence of flammable anesthetic mixture with air, oxygen or nitrous oxide.
- Follow local ordinances or regulations for proper disposal of device, accessories and packaging.
- Caution: Federal law requires this device to be sold by a physician or by the order of a physician.


WARNING

Read the following treatment warning protocol instructions carefully before operating the cold rush cold therapy system.

- The flow control knob should never be set at the maximum settings (number “4” or “MAX”) for more than 30 minutes.


or

- The Cold Rush pad should never be applied directly to the patient's skin. Do not use the Cold Rush Cold Therapy System without first applying an insulation barrier between the Cold Rush pad and the

patient's skin.

- Patients should never utilize cold therapy while they are asleep or unable to check their skin every 1-2 hours unless specifically instructed to do so by the prescribing health care provider.
- Patients must check their skin every 1-2 hours once steady flow of water through the pad is achieved. Turn the flow control knob counterclockwise to increase the temperature (warmer), and clockwise to decrease the temperature (colder).
- When turning the unit on, turn the flow control knob clockwise to the number "3" setting wait 10 minutes to allow the temperature to stabilize before making any additional temperature adjustments.

- The following patients should ALWAYS turn the flow control knob counterclockwise to the number "7" setting when using cold therapy
 - Patients whose treatment area has been desensitized due to local anesthesia or regional nerve blockers
 - Patients taking prescribed medications that may negatively affect peripheral vascular circulation, including beta adrenergic blockers and local epinephrine use (such as in local anesthetics)
 - Patients taking prescribed medications that negatively affect their mental capacity or judgment
 - Patients suffering from diabetes
 - Patients who have recently undergone hand/wrist or foot/ankle surgery
 - Patients suffering from cognitive disabilities
 - Young children and the elderly
- Caution: Contact your healthcare provider immediately and discontinue the use of your Cold Rush Cold Therapy Unit if you experience any adverse reaction to the treatment area such as: increase in pain, burning, blistering, itching, increased swelling, discoloration of the skin, increased redness, welts or any other changes in the appearance of the skin.
- Make sure that the insulation barrier is dry.
- Care should be taken during prolonged use for children, diabetics, incapacitated patients and those with decreased skin sensitivity or poor circulation.
- A licensed healthcare practitioner must determine the correct frequency and duration of treatment for each patient.
- If prescribed by a licensed healthcare practitioner for home use, the healthcare practitioner must provide each patient with training for appropriate use and proper application of the Cold Rush Cold Therapy System.
- The health care practitioner must monitor the patient's use of the device to assure appropriate use and proper application of cold therapy, including but not limited to application of an insulation barrier between the patient's skin and the Cold Rush Pad.
- The product should never be left unattended when plugged in.
- Never drop or insert any object into any opening or hose.
- Do not operate where aerosol (spray) products are being used or where oxygen is being administered.
- Do not attempt to drag or reposition the Cold Rush Cold Therapy unit

by pulling or tugging on the insulated hoses, as this may cause damage to the unit and/or cause the hoses to disconnect. Use handle when moving the Cold Rush Cold Therapy unit.

- If you experience disturbance in any other electrical equipment near the Cold Rush unit, please try to unplug it from the wall socket to verify that it is not the source of disturbance.

ELECTRIC SHOCK HAZARD

- Unplug unit before filling it with ice and water. To avoid danger of electric shock be certain that hands are dry before inserting or removing the power supply from wall socket.
- Avoid getting water or ice on the flow control panel or power adapter.
- Always unplug this product immediately after each use.
- Do not use while bathing.
- Do not place or store product where it can fall or be pulled into a tub or sink.
- Do not place in or drop into water or other liquid.
- Do not reach for a product that has fallen into water. Unplug immediately.
- Never operate this product if it has a damaged electrical cord or plug, if it is not working properly, if it has been dropped or damaged, or dropped into water. Return the product to your healthcare provider for examination and repair.
- Keep cord away from heated surfaces.

NOTE: To create the safest operating conditions possible, Össur has designed the Cold Rush Cold Therapy System to operate on a 12 volt power supply, significantly reducing the chance of dangerous electric shock. Great care should always be taken when using electricity near water.

PRODUCT DEFINITIONS

- The Cold Rush Cold Therapy System is designed for single patient use only.
- Use of Cold Rush Cold Therapy System is designed to reduce localized pain and swelling following surgical procedures or trauma.

LIST OF PARTS (Figure 1)

- a. Cold Rush Cold Therapy System
- b. Hose
- c. Couplings
- d. Flow Control Knob
- e. DC Connector
- f. 12V Power Supply
- g. Cold Rush Pad (Sold Separately)

TECHNICAL SPECIFICATION

Do not use any other power supply other than the one provided with the unit.

ELECTRICAL - POWER SUPPLY

- Model no: UES06WU-120050SPA
- Input: 120V US, 240V EU
- Input frequency: 60Hz US, 50Hz EU
- Output power: 6W
- Output Voltage: 12V

- Output Current: 0.5A

SAFETY STANDARDS

Cold Rush™ Cold Therapy System is tested and certified to comply with the IEC60601-1 standard of electrical safety of medical devices and IEC/EN60601-1-2, electromagnetic compatibility for medical electrical devices. The company fulfills the requirements of ISO13485.

SETUP AND OPERATION

1. Ensure cold therapy unit is turned off, turn flow control knob counter clockwise to the off position. Do not plug in to wall socket **(Figure 2)**
2. Fill reservoir with water to the indicated fill line **(Figure 3)**
3. Add chopped or cubed ice to the reservoir up to the indicated fill line **(Figure 4)**
4. Attach the provided 12V power supply to DC power port, do not plug into wall circuit **(Figure 5)**
5. Apply insulation barrier over the patient's skin in the area to be treated **(Figure 6a)**
6. Apply Cold Rush Pad over insulation barrier. The Cold Rush Pad should never be applied directly over the patient's skin **(Figure 6b)**
7. Connect Cold Rush pad to hose extending from the unit. To ensure proper connection, push couplings together until they snap into place **(Figure 7)**
8. Plug the provided 12v power supply to wall socket (120V US, 240V Europe) **(Figure 8)**
9. Turn flow control knob clockwise to the number "3" setting, this will turn unit on **(Figure 9)**
10. Wait 10 minutes to allow temperature to stabilize before adjusting the temperature for comfort **(Figure 10)**
 - Turn flow control knob counterclockwise to increase temperature.
 - Turn flow control knob clockwise to decrease temperature.
11. To turn the cold therapy unit off, turn flow control knob counterclockwise until it stops.
12. Always turn the cold therapy unit off before disconnecting the Cold Rush pad or any hoses. To disconnect the hoses, press the release tabs on the coupling and pull the connectors apart.
13. Drain all of the water out of the cold therapy unit after each use.
14. Follow directions 2-10 before each use.

MAINTENANCE

Do not try to repair your Cold Rush unit as this may create a user hazard and will automatically void warranty. After talking with your healthcare provider, to ensure you are following the operating instructions properly, the unit shall be sent to Össur for repairs or maintenance. No other company or individual is allowed to service or repair the unit in any way.

CLEANING AND PROLONGED STORAGE

- When you have finished all of your cold therapy treatments, prior to prolonged storage, fill the unit with one gallon of water (or enough so the water level is at the ice fill line.
- Add 2 tablespoons of liquid bleach or 3 ounces of 3% hydrogenated peroxide to the water and mix.
- Connect hoses to the Cold Rush Pad, turn unit on to MAX setting.
- Allow water mixture to circulate for 5 minutes.
- Turn unit off and unplug the power supply.
- Carefully drain the water mixture out of the unit and store in a cool

dark place leaving the lid slightly open. (Exposure to sunlight and extreme heat may damage hoses and unit.)

- When reusing the unit after storage, fill the reservoir with clean water and turn the pump on to circulate the water through the hoses and flush the remaining water mixture out of the system (use old pad for this procedure, discard the pad after cleaning is completed).
- Dispose water and refill the reservoir with fresh ice and water combination to begin cold therapy treatment. Follow operating instructions 2-10 before each use.
- Before you use the unit, inspect hoses and new pad for any tears or breaks. Contact your healthcare provider to obtain new Cold Rush pad.
- Clean or replace insulation barrier as recommended by your healthcare practitioner.

ENVIRONMENTAL

- Operating Temperature: +10°C to 40°C
- Operating Humidity: 30% - 75% Relative Humidity
- Atmosphere Pressure: 700h Pa – 1060 h Pa
- Shipping and Storage Temperature: -20°C to 60°C
- Shipping and Storage Humidity: 10% - 90% Relative Humidity, non-condensing

ACCESSORIES

The following accessories are available for your Cold Rush Cold Therapy System at an additional charge. To order, contact your healthcare provider.

- Cold Rush Ankle Pad
- Cold Rush Hip Pad Left
- Cold Rush Hip Pad Right
- Cold Rush Knee Pad
- Cold Rush Lumbar Pad
- Cold Rush Rectangular Pad Small
- Cold Rush Shoulder Pad Small
- Cold Rush Shoulder Pad Large
- Cold Rush Universal Pad

LIMITED WARRANTY

Össur warrants, to the original purchaser, that this product is free from defects in materials and workmanship for six: 6 months from initial purchase. This warranty does not apply if this product: 1 was not purchased from Össur or its authorized dealer; 2 has been altered in any way, or: 3 has not been used in accordance with Össur's Instructions For Use (IFU). This Limited Warranty does not cover damage due to accidents, neglect, misuse, or operation beyond capacity, parts damaged by improper installation, substitution of parts not approved by Össur, or any other alteration or repair by others that, in Össur's judgment, materially or adversely affects the product or part.


For complete warranty guidelines please visit www.ossur.com

ELECTROMAGNETIC COMPATIBILITY PRECAUTIONS AND WARNINGS

The Cold Rush Therapy System needs special precautions regarding electromagnetic compatibility (EMC). Specifically it needs to be installed and put into service according to the EMC information provided as follows:

- The Cold Rush Therapy System should not be used adjacent to or stacked with other equipment. In case adjacent or stacked use is necessary, the Cold Rush Therapy System should be observed to verify normal operation in the configuration in which it will be used.
- The Cold Rush Therapy System may be susceptible to electromagnetic interference from portable and mobile RF communications devices such as mobile (cellular) telephones.
- The Cold Rush Therapy System may be interfered with by other equipment, even if that other equipment complies with CISPR EMISSION requirements.

Guidance and Manufacturer’s Declaration – Electromagnetic Emissions		
The device is intended for use in the electromagnetic environment specified below. The customer or user of device should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic Environment – Guidance
RF Emissions CISPR 11	Group 1	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
RF Emissions CISPR 11	Class B	
Harmonics Emissions IEC 61000-3-2	Class B	
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration – Electromagnetic Emissions			
The device is intended for use in the electromagnetic environment specified below. The customer or user of device should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete, or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, Short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5 % U_T (>95 % dip in U_T) for 0.5 cycle. 40 % U_T (60 % dip in U_T) for 5 cycles. 70 % U_T (30 % dip in U_T) for 25 cycles. <5 % U_T (>95% dip in U_T) for 5 s.	<5 % U_T (>95 % dip in U_T) for 0.5 cycle. 40 % U_T (60 % dip in U_T) for 5 cycles. 70 % U_T (30 % dip in U_T) for 25 cycles. <5 % U_T (>95% dip in U_T) for 5 s.	Mains power quality should be that of a typical commercial or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U_T is the a.c. mains voltage prior to application of the test level.			

Recommended separation distances between portable and mobile RF communications equipment and the device.			
The device is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the device can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz d=1,2√P	80 MHz to 800 MHz d=1,2√P	800 MHz to 2,5 GHz d=2,3√P
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.			
NOTE 1	At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.		
NOTE 2	These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		

Guidance and Manufacturer's Declaration – Electromagnetic Emissions			
The device is intended for use in the electromagnetic environment specified below. The customer or user of the device should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment – Guidance
			Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance:
Conducted RF IEC 61000-4-6	3 Vrms, 150 kHz to 80 MHz	3 Vrms	$d=1.2\sqrt{P}$
Radiated RF IEC 61000-4-3	3 V/m, 80 MHz to 2.5 GHz	3 V/m	$d=1.2\sqrt{P}$ 80 MHz to 800 MHz $d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz
			where “P” is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey*, should be less than the compliance level in each frequency range**. Interference may occur in the vicinity of equipment marked with the following symbol: 
Note 1	At 80 MHz and 800 MHz, the higher frequency range applies.		
Note 2	These guidelines may not apply in situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.		
*	Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.		
**	Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.		

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