

Specifications	Box	Assembled Product
Weight	12.1 lbs	8.0 lbs
Length	19.7 in	18.1 in
Width	19.2 in	17.0 in
Height	5.1 in	11.4 in

Limited Lifetime Warranty

Your Drive Medical brand product is warranted to be free of defects in materials and workmanship for the lifetime of the original consumer purchaser. This device was built to exacting standards and carefully inspected prior to shipment. This Lifetime Limited Warranty is an expression of our confidence in the materials and workmanship of our products and our assurance to the consumer of years of dependable service.

In the event of a defect covered by this warranty, we will, at our option, repair or replace the device.

This warranty does not cover device failure due to owner misuse or negligence, or normal wear and tear. The warranty does not extend to non-durable components, such as rubber accessories, casters, and grips, which are subject to normal wear and need periodic replacement.

If you have a question about your Drive device or this warranty, please contact an authorized Drive dealer.



© 2015 Medical Depot, Inc. All rights reserved.
Drive is a trademark of Medical Depot, Inc.
Port Washington N.Y. 11050 USA Made in China

Oxygen Cylinder Metal Stand

item #:
18149KD



www.drivemedical.com



Product Identification

For use with M60, M90, M120, H and T Cylinders from 7" to 9" in diameter.

Please be sure all items are included

1. One stand body
2. One adjustable reversible bracket
3. Three hand knobs



Product Assembly

CAUTION Proper assembly required to avoid tipping hazard. Be sure to follow assembly instructions carefully.

1. Remove the 18149KD from all of the packaging material and verify all pieces are included.
2. Pull the side post and horseshoe shaped assemble toward the front of the stand until the side posts are perpendicular to the floor.
3. Swing the horseshoe bracket so the bolt faces the rear of the stand and is parallel to the floor.
4. Move the rear support bar into position and insert the bolt from the horseshoe shaped assembly into the hole in the rear support. Add and tighten the hand knob.



5. Insert the cylinder until the rear is flush with the back of the horseshoe bracket.
6. Attach the front adjustable reversible bracket, according the size of the cylinder, and tighten each hand knob.



Larger Cylinders



Smaller Cylinders