

Be Independent.™

# 30" SAFETY BED RAIL™ & PADDED POUCH

BECAUSE YOUR BED RAIL CAN PIVOT OUT OF THE WAY



PRODUCT #8051



WATCH A DEMO





Getting Out of Bed is Easier With the Hand Rail



Includes a Bonus Organizer Pouch



Safety Bed Rail Pivots Out of Way



Safety Strap Holds You and the Rail in Place



## FEATURES & BENEFITS

- **Bed Rail/Hand Rail** - Prevents falling from bed and can be used for assistance to sit up in bed.
- **Padded Pouch** - Safety bumper pad with 4 convenient pockets per side.
- **Collapsible** - Rail pivots down to side of bed to allow user to get out of bed or provide space for making bed.
- **Reversible** - Can be placed on either side of the bed.
- **Dual Safety Strap** - Secures bed rail to any sized home or hospital bed frame.

## SPECIFICATIONS

*Width of Rail: 30"*

*Height of Rail: 23" high from the base of mattress*

*Gap Between Rails: 3"*

*Bed Type: Attaches to any home or hospital bed with included safety strap*

*Weight of Product: 12 lbs*

*Weight Capacity: 400 lbs*

*Assembly: Installs in minutes with 4 bolts and allen wrench (included)*

*Package Dimensions: 32" L x 21.5" W x 1.25" H 12 lbs*

*Case Dimensions: 32" L x 21.5" W x 5"H; 50 lbs*

## EASY ASSEMBLY (No Additional Tools Required)



Attach L Shaped Tubes



Slide on Strap



Attach Rear Crossbar and Padded Pouch



Secure to Bed Frame with Included Safety Strap