



# Respiratory Care

Breathe Easy...A Complete  
Line of Respiratory  
Products





# Ensuring Your Success in Respiratory and Anesthesia Care

Medline is a market leader in  
one-source respiratory care for

Oxygen Therapy, Suction, Trach Care,  
Asthma Management, CPAP and Aerosol  
Therapy. We are clearly focused on the  
future of Respiratory products and strongly  
believe service and quality are the keys to  
daily success in this business.

Our focus will be to listen  
to you about your needs and  
partner with you to be your  
one BEST source of quality  
products, service and  
programs.

## Table of Contents

<b>Oxygen Therapy</b>	<b>3</b>
Oxygen Concentrators	
AMBU Bags	
Mouth-to-Mouth Resuscitator/Face Shields	
Oxygen Cannulas & Tubing	
Oxygen Masks	
Oxygen Connectors and Adaptors	
Pulse Oximeters & Accessories	
Oxygen Analyzer	
Oxygen Regulators and Conservers	
Oxygen Flow Meters	
Oxygen Cylinders and Accessories	
<b>Aerosol Therapy</b>	<b>17</b>
Large Volume Air Compressor	
Large Volume Nebulizer	
Corrugated Tubing & Drainage Bags	
Aerosol & Trach Masks	
<b>Asthma Management</b>	<b>20</b>
Nebulizer – Compressors	
Peak Flow Meters	
Holding Chambers	
Incentive Spirometers	
Nebulizer Kits	
Saline Vials	
<b>Suction</b>	<b>26</b>
Suction Pumps	
Open Suction Catheters / Kits	
Closed Suction Systems	
Yankauers, Poole Tips & Mucus Traps	
Suction Canisters & Absorbent Products	
Suction Tubing	
Suction Regulators	
Suction Irrigation Solutions	
Solidifiers	
<b>Tracheostomy Care</b>	<b>40</b>
Care & Cleaning Trays	
Trach Tube Holders	
Trach-Vox	
<b>CPAP/BiPAP Therapy</b>	<b>43</b>
CPAP/BiPAP Machines	
Masks & Nasal Pillow	
Headgear	
Tubing	
Mask Cleaner	
<b>Ventilation</b>	<b>46</b>
Transport Ventilation	
Omni-Vent	
Auto-Vent	
i-Vent	
Heat Moisture Exchangers ( HME )	





## Oxygen Therapy

Patient comfort, clinical and customer effectiveness inspire our products. Medline's extensive line of Oxygen Therapy products increase the amount of oxygen in the blood and help improve patient comfort. Our products range from our exclusive oxygen concentrator to our wide variety of O<sub>2</sub> delivery regulators, conservers, humidifiers, disposables and a high-quality AMBU bag line. We also offer the unique Soft-Touch cannula, that is available in a variety of sizes and styles. It is a soft, pliable, over-the-ear nasal cannula that provides a better patient fit and comfort.

# Compact Design & Solid Engineering converge in the NEW EndurO<sub>2</sub> Oxygen Concentrator

## Design keeps portability in mind

- Compact, continuous 5-liter concentrator
- Weighs only 41 lbs.
- Wide-set wheels for stability
- Integrated handles on top and on each side


## Maintenance is practically effortless

- Three-year limited warranty
- Easy access to all three filters without tools
- One way check valve prevents backflow and allows every 4 year bacterial filter change
- Easy access to all internal components
- Nationwide warranty/service center system
- Hotline: 1-866-225-9804
  - maintenance/repair
  - technical questions
- Product training available
  - \*basic operations and maintenance DVD
  - in-depth routine service and maintenance

## Convenience is built right in

- Easy view instrument panel (65° angle)
  - power switch
  - reset button (circuit breaker)
  - flow settings
  - hour meter
  - indicator lights
- Filter replacement frequency
  - cabinet (one/two times a year)
  - compressor (annually)
  - bacterial/viral (every four years)

**\*Basic operation and maintenance video is easy to understand.**



**Available from your Medline sales rep.**



Specification	Model HCS021	Model HCS02
Equipped with OCI (Oxygen Concentrator Indicator)	No	Yes
Alarms	Power Failure High and Low Pressure	Power Failure High and Low Pressure Low O <sub>2</sub> Purity
Dimensions	22"H x 15"W x 12"D	
Weight	41 lbs. (18 kg)	
Electrical	115 V AC, 60 Hz, 3.8 Amps, 420 W max. Double insulated cabinet	
Max. Outlet Pressure	8.5 +/- 0.5 psi (0.6 +/- 0.1 bar)	
Flow Range	1 - 5 L.P.M.	
Oxygen Concentration	93% +/- 3%	
Filter	Cabinet Compressor Intake Bacteria	
Sound Level	50 dBA +/- 2 dBA	
Operating System	Time Cycle / Pressure Swing	
Warranty	Three-year limited warranty	
Operating Environment	Ambient Temperature 50°F to 95°F (10°C to 35°C) Atmospheric Pressure: Sea Level to 7000 ft.	
Storage/Transportation	Temperature: -40°F to 158°F (-40°C to 70°C) Humidity: 10% to 100% (including condensation) Atmospheric Pressure: Sea Level to 7000 ft.	

Item #	Description	Pkg
HCS021	Oxygen concentrator without low-O <sub>2</sub> monitor	each
HCS02	Oxygen concentrator with low-O <sub>2</sub> monitor	each
AGICABINET	Replacement cabinet filter	10/pk
AGIINTAKE	Replacement intake filter	10/pk
AGIBACTERIA	Replacement bacterial filter	3/pk



## Oxygen Analyzers

Max-O<sub>2</sub> features an integrated O<sub>2</sub> Sensor for greater reliability, and is ideal for homecare and nursing home use. Includes corrugated tubing adapter.

The “**Handi**” is designed to be maintenance-free and disposable at the end of the sensor life. This analyzer is great for spot checking O<sub>2</sub> percentage in hospitals, clinics, and in home care environments.

Both feature

- Easy, one-step calibration
- 2-year warranty—analyzer and sensor

Item #	Description	Pkg
<b>HCSR217P62</b>	Max-O <sub>2</sub>	each
HCSR217P32	Carry bag for Max-O <sub>2</sub>	each
HCSR217P11002	Back Assembly	each
HCSR2117P23	Dovetail Assembly	each
HCSR206P75	Pole Clamp	each
<b>CERR207P01</b>	Handi Mini O <sub>2</sub> Analyzer	each



Max-O<sub>2</sub> Oxygen Analyzer  
HCSR217P62



CERR207P01



HCS64184

## Mouth to Mask Resuscitator Protects the Rescuer and Victim with a Highly Effective Viral and Bacterial Filter

This ready-to-use CPR mask is made with protection in mind as it features a one-way filtering system to protect the victim and rescuer at the same time. The clear plastic design allows for easy placement while the soft inflated cuff provides patient comfort and a better fit. Standard 22mm fittings allow the mask to be used with any resuscitator. Simplified graphic and anatomical illustrations are printed on the hard plastic carrying case.

Item #	Description	Pkg
HCS64184	Mouth to Mask Resuscitator w/ Valve	12/cs

## CPR Barrier Mask Is a Must-Have for All Health Care Workers

This small and flexible face shield features a one-way valve and filter to help protect the rescuer and victim alike. Simple graphic and anatomical illustrations are printed on the face shield to help ensure proper placement and use.

Also available with a key ring pouch for immediate access in case of an emergency.

Item #	Description	Pkg
HCS63002	CPR Barrier Shield	50/cs
AMB248201103	CPR Barrier Shield with key ring pouch	each



HCS63002



AMB248201103



## Ambu SPUR II Manual Resuscitators

Ambu SPUR II features a unique shape that is small and easy to hold, yet provides optimum stroke volumes with just one hand. The bag itself is thin and responsive and the SafeGrip surface ensures that the Ambu SPUR II can be firmly squeezed.

All Ambu SPUR II resuscitators come in individual, re-sealable plastic carrying bags, complete with one or more masks and many special accessories. The bags themselves are color-coded for fast identification of infant, pediatric, and adult sizes.

### KEY FEATURES

- Unique single-shutter valve system for reliable functionality
- Integrated handle for user comfort and uniform compression
- Thin-walled compression bag allows for lung compliance and “feel”
- SafeGrip surface for secure handling in stressful environments
- Medication port permits quick medication delivery without disconnecting the Ambu SPUR II from the ET tube

Item #	Description	Pkg
<b>AMB520611000</b>	Adult Bag w/Med Adult Mask and Bag Reservoir	6/cs
<b>AMB530613000</b>	Pediatric Bag w/Toddler/Ped Mask and Bag Reservoir	6/cs
<b>AMB544212000</b>	Infant Bag w/Infant Mask and Tube Reservoir	12/cs
AMB520211000	Adult Bag w/Med Mask, Syringe Port and Bag Reservoir	12/cs
AMB521611000	Adult Bag w/Med Mask, Syringe Port, Bag Reservoir and 40" Tubing	6/cs
AMB530613000	Pediatric Mask w/Toddler/Ped Mask and Bag Reservoir	6/cs
AMB531613000	Pediatric Bag w/Toddler/Ped Mask and 40" Tube Reservoir	6/cs
MSN540212000	Infant Bag w/Infant Mask and Bag Reservoir	12/cs

### Adult Specifications

Stroke volume: One hand	800ml
Two hands	1100ml
Resuscitator volume:	1475 ml
Dimensions (length x diameter):	295 x 127 mm
Weight (Reservoir O2 tube and mask):	345 g (40") 397 g (60")
Dead space:	<6 ml
Tube reservoir: 40"	
Diameter	22 mm (7/8")
Volume	393 ml
Optional: 60"	
Diameter	22 mm (7/8")
Volume	590 ml

### Specifications All Versions

Patient connector:	
Outside:	22 mm male
Inside:	15mm female
Expiratory connector (for PEEP valve attachment):	30 mm male
Forward and backward leak:	Not measurable
MediPort:	Standard Luer LS 6
Recommended operating temperature:	0 F to +122 F (-18 C to +50 C)
Storage:	Tested at -40 F and +140 F (-40 C and + 60 C)

## Medline Nasal Cannula — Comfortable, Economical and Latex-Free

Medline nasal cannulae are latex-free and available in numerous sizes, shapes and lengths. Our Soft-Touch cannulae are second to none in patient comfort when caring for long term oxygen patients.



HCS4514 (Soft-Touch)



\*HCS4511

Item #	Description	Pkg
<b>HCS4514</b>	<b>Adult Cannula, Soft-Touch, 7' crush resistant tubing</b>	<b>50/cs</b>
HCS4510	Adult Cannula Only, Soft-Touch	50/cs
<b>HCS4511</b>	<b>Adult Cannula, Economy, 7' crush resistant tubing*</b>	<b>50/cs</b>
HCS4504	Adult Cannula, Soft-Touch, 4' crush resistant tubing	50/cs
HCS4516A	Adult Cannula, Soft-Touch, 14' crush resistant tubing	50/cs
HCS4515A	Adult Cannula, Soft-Touch, 25' crush resistant tubing	25/cs
HCS4517	Infant Cannula, Soft-Touch, 7' crush resistant tubing	50/cs
HCS4518	Pediatric Cannula, Soft-Touch, 7' crush resistant tubing	50/cs
WST0101	Mask, Nasal, BI Flo Cannula	25/cs
HUD1850	Dual Lumen Cannula, 5' standard tubing	25/cs
WST0137	Adult Micro Cannula, 7' standard tubing	50/cs
HCS0440	Ear Cushion Set	50/cs

**All cannulae feature crush-resistant tubing as a standard feature for greater patient safety.**

**Our “Soft-Touch” Cannulas have curved, non-flared nasal prongs**

**\*Our economy cannula is ideal for the short term post-operative patient needs and even in some long term care cases. It has straight, non-flared nasal prongs**

## Medline Oxygen Tubing

Medline oxygen extension tubing comes standard, latex-free, crush resistant, and low memory. Our Unique “relaxed” low memory tubing helps prevent undo kinking. Violet relaxed tubing provides greater visual contrast for improved safety. Available in a variety of lengths.

Item #	Description	Pkg
HCS4507	7' crush resistant, low memory tubing, clear	50/cs
HCS4507V	7' crush resistant, low memory, violet color tubing	50/cs
HCS4524	14' crush resistant, low memory tubing, clear	50/cs
HCS4525	25' crush resistant, low memory tubing, clear	25/cs
HCS4525V	25' crush resistant, low memory, violet tubing	25/cs
HCS455010	50' crush resistant, low memory tubing, clear	10/cs
HCS4550V	50' crush resistant, low memory, violet tubing	20/cs



HCS4507V





## Medline Disposable Oxygen Masks Fits All Tubing

Latex-free and available in a variety of styles, each mask is designed for patient comfort with a soft, anatomical form. Barbed fittings help prevent patient from disconnecting the mask while the elastic straps and adjustable nose clips secure it for a better fit.

Item #	Description	Pkg
<b>Adult</b>		
HCS4600	Medium Concentration Mask, 7 ft. tubing	50/cs
HCS4610	High Concentration Mask, Reservoir Bag, 7 ft. tubing	50/cs
<b>HCS4640</b>	Non-Rebreather Mask, Reservoir Bag, Safety Vent, Check Valve, 7 ft. tubing	50/cs
<b>HSK120E</b>	Variable Select Mask, 6 Color-Coded Diluters from 24%–50%, humidity hood, 7 ft. tubing	50/cs

<b>Pediatric</b>		
HCS4601	Medium Concentration Mask, 7 ft. tubing	50/cs
<b>HCS4642</b>	Non-Rebreather Mask, Reservoir Bag, Safety Vent, Check Valve, 7 ft. tubing	50/cs
HSK227	Variable Select Mask, 6 Color-Coded Diluters from 24%–50%, humidity hood, 7 ft. tubing	50/cs



HCS4642



HCS4640



HSK120E


Detail of  
Humidity connector



Connectors & Adapters

Medline offers a variety of common and not so common connectors for different Respiratory and Anesthesia applications.

Item #	Description	Pkg
HCS0590	Connector, oxygen tubing, 360 degree swivel	25/bx
HCS65100	Connector, oxygen tubing	200/cs
HSK962A	Connector, 22mm OD and 15mm OD/22mm OD, for large bore tubing	50/cs
HCS1642	Adapter, Pressure line, also ideal for O2 entrainment	50/cs
WST219	Adapter, Aerosol Tee	50/cs
AGISH749	Connector, Oxygen tubing Splitter	5/pkg
HCS4541	Corrugated tubing, 6"	50/cs
HCS4555	Nut and Nipple Adapter (Christmas Tree)	50/cs
HCS2503	Disposable Water Trap	25/cs
HCS0700	Cylinder Wrench, small	50/cs
HCS166840	Brass Washers	50/bx



Medline Bubble Humidifiers Come in Several Sizes and Configurations

For respiratory therapy, with or without adapters. For use with 040 Humidifier adapters

Item #	Description	Pkg
HCS00340	Humidifier with Adapter, 340ml, Prefilled	20/cs
HCS00640	Humidifier with Adapter, 640ml, Prefilled	20/cs
HCS64377	Humidifier with Adapter, Empty, 6psi, Pop-off	50/cs



BCI3420Y

## Digit® Finger Oximeter

The Digit Finger Oximeter is a Pulse Oximeter and Sensor in one unit providing fast, reliable SpO<sub>2</sub>, pulse rate and pulse strength measurements on patients from pediatric to adult. Its pocket-size with a large LED display, one button keypad. The auto power shutdown after 8 seconds when not in use, makes it ideal for emergency, hospital or clinical environments, or for home use.

Item #	Description	Pkg
BCI3420Y	Digit Finger Oximeter, Yellow	each

**Great for Spot-Checking & Continuous!**

## Handheld Pulse Oximeter

Compact and lightweight, our handheld pulse oximeter is perfect for non-invasive spot-check measurement of oxygen saturation of arterial hemoglobin. It features automatic standby, power off, and record management. It can store up to 100 patient IDs and 200 records. Downloadable to computer. Multiple Sensor Capability.

Item #	Description	Pkg
<b>HCSPM50</b>	Pulse Oximeter, (Uses 4 AA batteries)	each
HCSPM50CM	Cable Data Link & Software	each
<b>HCSMIND50CSWC</b>	Pouch, protective w/ belt clip	each
<b>HCSR135P02</b>	Pulse Oximeter Ear Probe (Nellcor)	each
<b>PCL0603120</b>	Adult Reusable finger clip/probe (Nellcor)	each
PCL0623120	Pediatric Reusable finger clip/probe (Nellcor)	each
PCL1003120	Adult Reusable soft sensor (Nellcor)	each
PCL1023120	Pediatric Reusable soft sensor (Nellcor)	each



## Pulse Oximeters with Alarms

*Ideal for overnight sleep screening.*

Item #	Description	Pkg
<b>BCI3303</b>	Continuous Handheld Oximeter with rechargeable battery & alarms, fits neonate to adult	each
<b>BCI3402</b>	MiniCorr Handheld Oximeter with battery or AC power & audible and visual alarms, fits neonate to adult	each



BCI3303



BCI3402





## Disposable Oximeter Probes Offer a One Year Warranty

Our new line of replacement sensors are all latex-free and are applied with three different adhesives.


- Micro-foam – best for easy sensor application and re-positioning
- Medplast – best for its maximum adhesive for long-term uses in one location
- Transpore – best for transparency with re-positioning


These disposable sensors are available for use with BCI, Nellcor (Non-Oxi-Max 9 Pin), Nonin and Ohmeda pulse oximeters and come in a variety of sizes. **Other probes available, contact your Medline representative**

Item #	Description	Pkg
<b>Applied with Micro-foam </b>		
PCL0202121	Adult Disposable finger probe (Nellcor Monitor)	24/cs
PCL0202301	Adult Disposable finger probe (BCI Monitor)	24/cs
PCL0202122	Pediatric Disposable finger probe (Nellcor Monitor)	24/cs
PCL0202302	Pediatric Disposable finger probe (BCI Monitor)	24/cs
PCL0202123	Infant Disposable finger probe (Nellcor Monitor)	24/cs
PCL0202303	Infant Disposable finger probe (BCI Monitor)	24/cs
PCL0202124	Neonatal Disposable finger probe (Nellcor Monitor)	24/cs
PCL0202304	Neonatal Disposable finger probe (BCI Monitor)	24/cs





Micro-foam

<b>Applied with Medplast </b>		
PCL0342121	Adult Disposable finger probe (Nellcor Monitor)	24/cs
HCS22211P	Adult Disposable finger probe (BCI Monitor)	24/cs
HCS22112P	Pediatric Disposable finger probe (Nellcor Monitor)	24/cs
HCS22212P	Pediatric Disposable finger probe (BCI Monitor)	24/cs
HCS22115P	Infant Disposable finger probe (Nellcor Monitor)	24/cs
HCS22215P	Infant Disposable finger probe (BCI Monitor)	24/cs
PCL0342124	Neonatal Disposable finger probe (Nellcor Monitor)	24/cs
HCS22216P	Neonatal Disposable finger probe (BCI Monitor)	24/cs



Medplast

<b>Applied with Transpore </b>		
PCL0322101	Adult Disposable finger probe (Nellcor Monitor)	24/cs
PCL0322301	Adult Disposable finger probe (BCI Monitor)	24/cs
PCL0322122	Pediatric Disposable finger probe (Nellcor Monitor)	24/cs
PCL0322302	Pediatric Disposable finger probe (BCI Monitor)	24/cs



Transpore

### Reusable

PCL0603120	Adult Reusable finger clip/probe (Nellcor)	each
PCL0623120	Pediatric Reusable finger clip/probe (Nellcor)	each
PCL1003120	Adult Reusable soft sensor (Nellcor)	each
PCL1023120	Pediatric Reusable soft sensor (Nellcor)	each

### NONIN

PCL0202701	Adult Disposable finger probe (Nellcor Monitor)	24/cs
PCL0202702	Pediatric Disposable finger probe (Nellcor Monitor)	24/cs
PCL0202703	Infant Disposable finger probe (Nellcor Monitor)	24/cs



Reusable

Medline is a leader in the respiratory market and we are proud to offer you a complete line of oxygen-therapy supplies and equipment. Our innovative products are designed for clinical effectiveness as well as patient comfort. Medline stands behind its products through enforcing strict quality guidelines. As we look to the future, our goal is to provide you with superior quality products, and programs to strengthen your business with unsurpassed customer service.



## The EconO<sub>2</sub>mizer II “Lightweight”

Medline's new single-lumen pneumatic conserver

Medline introduces its new single-lumen pneumatic conserver – the EconO<sub>2</sub>mizer II Lightweight. This lightweight and compact unit has the greatest number of settings for standard continuous and conserving flows. The EconO<sub>2</sub>mizer II features all brass in high-pressure zones, which is recommended by the FDA for fire prevention.

- Small size and light weight make unit less obtrusive and more portable (L: 4.8", W: 2.4". H: 1.5", **Weight: 13 oz.**)
- Six conserving flow settings
- **Three continuous flow settings**
- Brass piston assembly – safer than aluminum
- Pulse design makes less noise and is less irritating to patient and those around them
- 5.7:1 conserving ratio throughout the entire range

Item Number	Description	Packaging
HCS2CON1L	Single-lumen pneumatic conserver, lightweight	each
HCS4514	Adult Cannula, Soft-touch, 7' crush-resistant tubing	50/cs
HCS4504	Adult Cannula, Soft-touch, 4' crush-resistant tubing	50/cs



# The EconO<sub>2</sub>mizer

Dual-lumen oxygen-conserving device

Medline's EconO<sub>2</sub>mizer represents a breakthrough in conserving-device technology. With components like a patented all-brass diaphragm, the EconO<sub>2</sub>mizer reduces your service calls and provides patients with a lighter and more compact design.

- Zero to 6 lpm and 11 flow settings
- Dual mode for conserving and continuous flow with 11 settings for each mode (No additional regulator required)
- 4:1 conserving ratio for fewer service calls
- Solid brass core for maximal safety
- Laser-etched markings prevent chipping and flaking



HCSO2CONS



Dual mode switch

Item Number	Description	Packaging
HCSO2CONS	EconO <sub>2</sub> mizer, conserver, 0-6 settings	each
HUD1850	Cannula, 5', dual lumen	25/cs

## Medline Electronic Mini-Conserver

The Mini-Conserver features a 5:1 conserving ratio and automatic shutoff. In continuous mode, it has a flow rate of 2 lpm. The Mini-Conserver has a six-month battery life (approximately 100-M6 cylinders) from one AA alkaline battery.

- High-pressure brass chamber
- Battery life indicator and battery test mode
- Measures 1.9" x 3.3" x 3.3" and weighs only 14 oz.!
- Easy-off handle
- 870-style post valve
- Conserve flow settings - 1, 2, 3, 4, 5



HCS1001R



Detail of flow settings




Detail of mode selector

Item Number	Description	Packaging
HCS1001R	Mini-Conserver	each
HCS4514	Adult cannula, soft-touch, 7' crush-resistant tubing	50/cs
HCS4504	Adult Cannula, Soft-touch, 4' crush-resistant tubing	50/cs

**NEW!**

## Low-Cost Regulators

- These regulators are fitted with a brass sleeve around key parts
- All regulators are laser-etched to eliminate flaking and ruboff.
- Test method for evaluating the ignition
- Sensitivity and fault tolerance of oxygen
- Regulators for all medical and emergency applications.



Item Number	Connection	Flow Range	Flow Increments
HCS8708M	870 CGA connection	0-8 lpm	0, .12, .5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 8
<b>HCS8715M</b>	<b>870 CGA connection</b>	<b>0-15 lpm</b>	<b>0, .12, .5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15</b>
HCS8725M	870 CGA connection	0-25 lpm	0, .5, 1, 2, 3, 4, 5, 6, 8, 10, 15, 25
<b>HCS5408M</b>	<b>540 CGA connection</b>	<b>0-8 lpm</b>	<b>0, .12, .5, 1, 1.5, 2, 2.5, 3, 4, 5, 6, 8</b>
HCS5415M	540 CGA connection	0-15 lpm	0, .12, .5, 1, 1.5, 2, 3, 4, 6, 8, 10, 15

## Medline Oxygen Regulators


All-brass high-pressure chamber means they're made to last

Medline oxygen regulators use an advanced, patented Intelliflo® disc to regulate oxygen flow down to 1/32 lpm with very tight tolerance.

- Brass core for maximum safety
- Permanent laser etchings on knob and body
- All-metal anodized aluminum construction
- Gauge protector to prevent damage

Item Number	Description	Packaging
<b>HCSB8708M</b>	<b>1/2 to 8 lpm, small tank regulator</b>	<b>each</b>
HCSB8715M	0 to 15 lpm, small tank regulator	each
HCSB8704	1/32 to 4 lpm, small tank regulator	each
HCSB5408	1/2 to 8 lpm, large tank regulator	each
<b>HCSB5415</b>	<b>0 to 15 lpm, large tank regulator</b>	<b>each</b>
HCSB5404	1/32 to 4 lpm, large tank regulator	each
HCS166840	* Brass seal washer for O <sub>2</sub> regulators	50/bx
*Recommended by NFPA		

Call for availability of DISS Outlet





## Oxygen Cart Systems

### Rack'n Roll Modular Cylinder

Makes the moving and storing of aluminum medical cylinders easier. This complete modular product easily meets your rack and cart needs.

- Heavy-duty rubber tires (removable for stationary cart)
- Telescoping, height-adjustable holders
- Easy-to-hold hand grip for extra comfort
- Each six cylinder cart can link to make larger system (12, 18, 24 cylinders)
- Adjustable holder height allows each tank tray to accommodate different sizes



HCSRNR2001W

**Handle  
adjusts from  
36"-42"!**

### D-E Cylinder Rolling Cart

E-Tank cart is for use with E and D-Tank gas cylinders, allowing a patient to easily tote oxygen.

Item Number	Description	Packaging
HCSRNR2001W	Rack'n Roll oxygen cart system	each
HCS53006	E-Tank oxygen cart with height adjustment	6/cs
HCS0700	E-Tank wrench (non-metal)	50/cs

HCS53006



D-Tank  
Height

E-Tank  
Height

## Oxygen Therapy Accessories

Buy all of your related components from one supplier!



LXF6M, LXF6C, LXF6D, LXF6E, LXF60

Item Number	Description	Packaging
HCSO2M6	Carry bag, shoulder style, for M6 cylinder	each
HCSO2C	Carry bag, shoulder style, for C cylinder	each
HCSO2D	Carry bag, shoulder style, for D cylinder	each
LXFM6★	Cylinder, oxygen, M6 w/o toggle, 870 CGA post valve	each
LXFMC★	Cylinder, oxygen, C w/o toggle, 870 CGA post valve	each
LXFMD★	Cylinder, oxygen, D w/o toggle, 870 CGA post valve	each
LXFME★	Cylinder, oxygen, E w/o toggle, 870 CGA post valve	each
LXFM60★	Cylinder, oxygen, 540 CGA post valve	each
LXFM22★	Cylinder, oxygen, M22 w/o toggle, 870 CGA post valve	each
HCSO2TOTE	6 pack carrier for M6 cylinders	25/cs

★Add T to end of item number if toggle valve is required.



HCSO2M6



HCSO2C



HCSO2D

**Great for M6 Cylinder  
Transporting!**



HCSO2TOTE



## Flowmeters available for both oxygen and medical air.

The durable polycarbonate hood and metering tube are impervious to the most common hospital cleaning and disinfecting solutions, and are rugged enough to stand up to the hospital/nursing home environment.

- All are 0-15 LPM flow range, with minimum 65 LPM flush flow rate.
- Most are available with either a brass or aluminum body (aluminum body flowmeters weigh almost half as much as the brass bodies)
- Calibrated at 50 psig and 70° F, and are 100% tested to ensure accuracy.
- Needle valves in the flowmeter bodies provide precise flow control.
- All flowmeters are back-pressure compensated.
- 5-year warranty

Weight (No adapter): Aluminum Body - 4.24 oz./122 gm  
Brass Body -6.64 oz./188 gm

Increments: 0 - 6 LPM in 1/2 LPM increments,  
6 - 15 LPM in 1 LPM increments

Flush Flow Rate: 65 LPM (depending on supply flow)

Max Inlet Pressure: 100 psig

PTO Flow Rate: Over 100 LPM (depending on supply flow)



GTFM197B15LCHP  
Oxygen Flowmeter with Power Take-off  
and Chemtron Adapter

GTFM187B15L00  
Air Flowmeter

OXYGEN FLOWMETERS, 0-15 LPM		
GTFM197A15L = Aluminum Body	Select aluminum or brass body	
GTFM197B15L = Brass Body		
GTFM197B15L00 = No Adapter	Select the adapter (connects to the wall outlet)	
GTFM197B15LOH = Ohio/Ohmeda Adapter		
GTFM197B15LDS = DISS Nut & Nipple Adapter		
GTFM197B15LCH = Chemetron Adapter		
GTFM197B-15LPB = Puritan-Bennett Adapter		
GTFM197B15LDH = DISS Hand-Tight Adapter		
GTFM197B15DHP = Power Take-Off (DISS male w/Check)	Add -P for power take-off assembly	

MEDICAL AIR FLOWMETERS, 0-15 LPM		
GTFM187A15L = Aluminum Body	Select aluminum or brass body	
GTFM187B15L = Brass Body		
GTFM187B15L-00 = No Adapter	Select the adapter (connects to the wall outlet)	
GTFM187B15L-OH = Ohio/Ohmeda Adapter		
GTFM187B15LDS = DISS Nut & Nipple Adapter		
GTFM187B15LCH = Chemetron Adapter		
GTFM187B15LPB = Puritan-Bennett Adapter		
GTFM187B15LDH = DISS Hand-Tight Adapter		
GTFM187B15DHP = Power Take-Off (DISS male w/Check)	Add -P for power take-off assembly	

Please contact your Medline sales representative or customer service for assistance with selecting your new flowmeters.

## Wall/Back Fittings




**CH**  
Chemetron® style fitting



**DH**  
DISS Hand-tight fitting



**DS**  
DISS Nut & Nipple fitting



**OH**  
Ohio (Ohmeda®)  
style fitting



**PB**  
Puritan-Bennett®  
style fitting





## Large Volume Aerosol Therapy

Medline's Aerosol Therapy product line covers a wide array of applications, from products that assist in care for the home tracheostomy patient, to a traditional post-operative patient who needs airway re-hydration. Medline can provide the equipment needed for any application. Whatever your aerosol needs, Medline has your solution.



## The 50 PSI MaxiCompressor

It is designed for continuous aerosol humidification applications and to run 24 hours a day, 7 days a week reliably. It features a locking pressure regulator valve, a removable protective lid, Thomas Compressor, and washable filter.

*The MaxiCompressor's excellent airflow characteristics mean lower operating temperatures and longer life.*

Item #	Description	Pkg
MIN501	MaxiCompressor	each

### \*Everything for the complete aerosol set-up



**NEW!**



HCS04128

## Medline Large Volume Nebulizers & Accessories Come in Several Sizes and Configurations to Fit Your Needs Precisely

For respiratory therapy, with or without adapters.

Item #	Description	Pkg
<b>HCS04128</b>	Large Volume Nebulizer w/Adapter, 1070ml, Prefilled Bottle	10/cs
HCS03128	Nebulizer cap for HCS04000	50/cs
HCS04000	1070ml Sterile Water; Bottle only	10/cs

## Corrugated Aerosol Tubing is Flexible and Accepts 22mm Adapters for a Secure Fit

For heated and non-heated use, with easy-fit tubing ends. Longer lengths are segmented every 6".

Item #	Description	Pkg
HCS4541	Corrugated Tubing, 6" length	50/cs
<b>HSK626</b>	Corrugated Tubing, 72" (6 ft.) length	50/cs
HSK600	Corrugated Tubing, 100 ft. Roll in Dispenser Box	each



HSK626

**Medline Aerosol Drainage Bags are Designed to Make Sure Patient Receives Mist Rather than Liquid**

Placed within the aerosol delivery system, the “T” or “Y” design traps water, ensuring that the patient does not receive liquid that has accumulated in the tubing.

Item #	Description	Pkg
<b>HSK974Y</b>	Aerosol Drainage System (Bag with Y adapter and hanger)	50/cs
HSK974	Aerosol Drainage System (Bag with T adapter and hanger)	50/cs
WST219	Aerosol Connector-Aerosol Tee	50/cs



**Tracheostomy Masks–Available in Adult and Pediatric Sizes**

For use in tracheostomy and laryngectomy aerosol therapy. Both sizes have swivel adapter and accept 22mm corrugated tubing.

Item #	Description	Pkg
<b>HSK245E</b>	Adult Mask with Swivel Adapter	50/cs
HCS4621	Pediatric Mask with Swivel Adapter	50/cs



**Medline Aerosol Masks—Manufactured to Provide the Ultimate in Patient Comfort**

Lightweight, latex-free masks are made of soft, comfortable vinyl.

Item #	Description	Pkg
<b>HCS4630</b>	Adult, Aerosol Mask, Elastic Strap	50/cs
HCS4631	Pediatric, Aerosol Mask, Elastic Strap	50/cs

<b>HCS4632</b>	Adult, Face Tent, Under Chin Design	50/cs
----------------	-------------------------------------	-------





# Asthma Management

The prevalence of asthma continues to rise in the U.S. and the demand for new innovative asthma products is rising as well. Medline takes great pride in offering you the most technologically-advanced asthma management products to care for your patients.

## Peak Flow Meters Include Daily Asthma Management Record to Help Patients Chart Their Progress

Built-in asthma management zone system offers a full range for both Pediatric or Adult use.

Item #	Description	Pkg
HCSE50104	Peak Flow Meter, Adult, 50-850 lpm	24/cs
HCSE50106	Peak Flow Meter, Adult, 50-850 lpm	96/cs
HCSE50204	Peak Flow Meter, Child, 50-350 lpm, 24/cs assorted colors	
HCSE50206	Peak Flow Meter, Child, 50-350 lpm, 96/cs assorted colors	
MLT3002052	Peak Flow Meter	10/cs
HDT755	Personal Best Peak Flow Meter	12/cs



**Great for kids and seniors!**



**NEW!**

## Leverhaler Features Easy One Hand Operation

A holding chamber with lever actuator to help users increase leverage on their MDI medications.

Item #	Description	Pkg
BSGLH1000	Leverhaler	10/cs

*Leverhaler pushes medication out for you*

## Opti-Chamber Provides Greater Drug Deposition

Assists in the delivery of Metered Dose Inhalers (MDI) and is ideal for geriatric and pediatric populations who have a difficult time coordinating inhaling MDI drugs. Can be used with any standard mask with 22 mm ID respiratory fitting, and comes standard with inspiratory flow signal.

Item #	Description	Pkg
HDT80010010	Opti-Chamber	10/cs
HDT80110	Opti-Chamber w/Mask, Small Newborn-18 months	10/cs
HDT80210	Opti-Chamber w/Mask, Medium 1 yr-6 yrs	10/cs
HDT80310	Opti-Chamber w/Mask, Large 6 yrs-Adult	10/cs







**NEW!**

## Medline's Aeromist™ Mini Nebulizer/Compressor System is Perfect for Home, Active Adult and Pediatric Patients

- Pressure:  $\geq 35$  psi
- Maximum Flow:  $\geq 13.5$  lpm maximum flow
- Treatment Time: 7-10 minutes
- Particle Size:  $< 5$  microns
- Total Weight: Only 3.75 lbs
- Dimensions: 10.4"L x 4.9"W x 3.7"D
- Frequency: 110V; 60Hz
- Grounding: Double-insulated cord and plug
- Integral nebulizer cup holder
- Easy access cap for inlet filter
- Short treatment times
- 5-Year warranty

**Comes complete with nebulizer kit, filters and accessory bag**

Item #	Description	Pkg
HCS50004	Aeromist Mini Plus Nebulizer/Compressor	6/cs
HCS50004RD	Aeromist Mini Plus Nebulizer/Compressor w/*Reusable and Disposable Neb Cup	6/cs



\*Reusable Neb Cup

***Private labeled nebulizers available  
Ask your Medline account rep about details.***

## Incentive Spirometer Helps Patients Perform and Monitor Their Own Post-Surgical Breathing Exercises

Medline's incentive spirometer makes it easy for patients and staff alike. One-way valve ensures patients inhale rather than exhale into the unit. Highly visible pistons and universal graphics (indicating correct inspiratory flow rate). Other features include an oxygen connector for supplemental oxygen, convenient handle, flexible popple tubing and bedrail holder.

Item #	Description	Pkg
DHD224000	4,000 ml Capacity and One-Way Valve	12/cs
DHD222500	2,500 ml Capacity and One-Way Valve	12/cs
DHD222500H	2,500 ml Capacity and One-Way Valve	each
DHD221200	Clini-Flow Low-Flow Spirometer	12/cs



DHD224000

DHD221200



**Dishwasher Safe!**



**Kido's Bear Provides Effective and Immediate Treatment for Pediatric Patients**

The Bear is a new innovative delivery system for aerosol and oxygen therapy for children with asthma or other acute respiratory diseases.

- Delivers medication from a standard nebulizer
- Provides oxygen (therapy only)
- Calms the pediatric patient while providing effective and immediate therapy
- A nebulizer and oxygen line are included with each bear
- Available in yellow and blue (three of each individually wrapped in a case)
- As effective as wearing an aerosol mask

Item #	Description	Pkg
HCSKBN4300	Kiddos Bear	6/cs



**Ideal for Pediatric & Alzheimer's Patients**

**TeddyNeb Nebulizer**

TeddyNeb is a pediatric nebulizer compressor. The bear design is less intimidating to children. Comes with a mouthpiece, 2 sizes of masks, 3 spare filters, and a back pack carry bag.

- Five year warranty
- Quiet performance
- Quick treatment time
- Safety shut off via an overload protector

Item #	Description	Pkg
HCSTB006	Teddy Nebulizer	6/cs
HCSTBFILTER	Spare Filters	5/pk



**New Design Coming Soon!**



## Aeromist™ Epic Nebulizer/Compressor System for Consistent Performance and Use in Any Home

Aeromist Epic is designed with a durable, piston-style compressor and provides treatments at 30-35 psi and 9 lpm. It features a thermal protector which helps prevent overheating. Weighing only 4.5 lbs., with electrical requirements of 110 volts and 60 Hz, it can be used safely just about anywhere. 3-year warranty.

Item #	Description	Pkg
HCS16004	Aeromist Epic Nebulizer/Compressor	6/cs

Ideal for those on the go!



## Sportmist Compressor Nebulizer System

Medline's portable Sportmist has everything needed for treatment-on-the-go, in an incredibly light package. The entire unit weighs only 20 oz. The Sportmist provides quick, reliable treatments with the optimal particulate size of 3-4 microns. All the items required for portability are included in the Nebulizer package. Each Sportmist comes with the battery, car adapter, AC adapter, carrying case, Nebulizer circuit, 2 masks (adult and pediatric), and replacement filters. The Sportmist provides a safe, convenient, and easy way to take treatments anywhere.

Item #	Description	Pkg
HCSKN9320	Portable Nebulizer Compressor with accessories	3/cs
HCSKN9320NOBT	Portable Nebulizer Compressor, no battery	each
HCSKN933001	Rechargeable battery	each
HCSKN932002	Replacement filters 5 inlet, 5 outlet	10/bx
HCSKN932004	Power adapter replacement	each
HCSKN932005	Adult mask with adapter	each
HCSKN932006	Pediatric mask with adapter	each
HCSKN932007	Car power adapter	each
HCSKN932008	Mouthpiece nebulizer	each
HCSKN932009	Drug ampoule assembly	each
HCSKN932010	Air tubing with adapter	each



Medline’s VixOne Nebulizer — One of the Best Performing Nebulizers on the Market — Offers a More Optimal Particle Size so the Patient Gets More Medicine

That means less residual medication, so no wonder it’s one of the best! Available in various kits to suit your patients’ needs.

Item #	Description	Pkg
<b>Disposable Nebulizer Kits</b>		
HCS4482	<b>Nebulizer, T Mouthpiece, 7 ft. Tubing</b>	<b>50/cs</b>
HCS4483	Nebulizer, T Mouthpiece, 7 ft. Tubing, 6" Reservoir	50/cs
HCS4485	Nebulizer, Adult Mask, 7 ft. Tubing	50/cs
HCS4486	Nebulizer, Pediatric Mask, 7 ft. Tubing	50/cs
HCS4492	Nebulizer, T Mouthpiece, 7 ft. Tubing, Side Stream	50/cs
HCS4493	Nebulizer, T Mouthpiece, 7 ft. Tubing, 6" Reservoir, Side Stream	50/cs
HCS4495	Nebulizer, Adult Mask, 7 ft. Tubing, Side Stream	50/cs
WST0312	<b>Nebulizer, Pediatric Animal Mask, 7 ft. Tubing</b>	<b>50/cs</b>
WST0382	<b>Nebulizer, Paci-Neb, 7 ft. Tubing</b>	<b>50/cs</b>

Reusable Nebulizer Kits

HCS301681	<b>Nebulizer, Reusable, mouthpiece, and 7ft. Tubing</b>	<b>20/cs</b>
HCS30168140	<b>Nebulizer, Reusable, adult mask, and 7ft. Tubing</b>	<b>20/cs</b>
HCS30168141	<b>Nebulizer, Reusable, peds mask, and 7ft. Tubing</b>	<b>20/cs</b>



HCS4482



WST0312



WST0382



HCS30168140



HCS30168141



HCS301681

Sterile Saline Solution Dispenses Quickly

Normal saline solution dispenses by the drop in 3, 5, and 15 ml increments, or quickly for procedural efficiency in tracheal cleansing.

Item #	Description	Pkg
HCS20039	Saline 0.9%, 3 ml, Unit Dose	1000/cs
HCS20039Z	Saline 0.9%, 3 ml, Unit Dose	100/bx
HCS20059	Saline 0.9%, 5 ml, Unit Dose	1000/cs
HCS20059Z	Saline 0.9%, 5 ml, Unit Dose	100/bx
HCS20069	Saline 0.9%, 15 ml, Unit Dose	144/cs



HCS20039





# Suction

Medline is a market value leader of suction products in the U.S. We manufacture, supply and directly support our full line of products that includes open and closed suction catheters, yankauers, suction canisters, absorbent pads, suction tubing, solutions and more. We also provide a full line of liquid medical waste solidifiers and disinfectants to minimize exposure to infectious waste.

# Medline Suction Catheters with the Contro-Vac™ Valve Helps Protect Caregivers

## Promote Better Infection Control

Why risk contamination when you can use Medline suction catheters with the exclusive Contro-Vac valve—proven more effective than chimney valves at preventing contact with fluids. Our 5-in-1 Connector and Stabilizing Bar ensure easy connection for any size tubing and a secure grip.



5-in-1 Connector & Stabilizing Bar

## Proven superior in lab tests

In an independent laboratory evaluation\*, Medline suction catheters with the Contro-Vac valve were compared to leading manufacturers' catheters with chimney valves. The results? Medline catheters were proven more effective at preventing the release of microbial contaminant during the suctioning process.

\* Geneva Lab Study (1995)



Medline's Contro-Vac valve




Traditional chimney-valve catheter

## A complete line of unique suction catheters and kits

Medline suction catheters are a better solution for your everyday needs. The following pages will reveal the value of our unique suctioning products.

**Medline  
Exclusive**



Easy to read  
depth markings


## Medline Depth-Marked Suction Catheters

The safe-suctioning technique, using depth-marked suction catheters, has been practiced for years with newborn and pediatric patients. The ATS (American Thoracic Society) strongly recommends the use of these depth-marked suction catheters for all patients\* to help reduce trauma and infection potential.

\*American Journal Respiratory Care Medicine  
Vol. 161, pgs 297-308

**Use our "Safe Suctioning Guide" as a new part of  
your Infection Control and Patient Safety Program**

**LIT340**







## Open Suction Catheter Kits (Dry)

- Full Range of Sizes (6-18Fr)
- Multiple Packaging Options

### Mini-Trays with Catheter & Gloves

- Smaller for Storage
- Tray is the Cup
- Best Price Option



Item #	Description	Pkg
DYND40986	6Fr (DeLee Tip), 2 Gloves	100/cs
DYND40988	8Fr (DeLee Tip), 2 Gloves	100/cs
DYND40980	10Fr (Whistle Tip), 2 Gloves	100/cs
DYND40981	12Fr (Whistle Tip), 2 Gloves	100/cs
<b>DYND48982</b>	<b>14Fr (Whistle Tip), Poly-Cath, 2 Gloves</b>	<b>100/cs</b>
<b>DYND40982</b>	<b>14Fr (Whistle Tip), 2 Gloves</b>	<b>100/cs</b>
DYND40983	16Fr (Whistle Tip), 2 Gloves	100/cs
DYND40984	18Fr (Whistle Tip), 2 Gloves	100/cs



DYND40982



DYND48982

### Sleeved Catheter & Cup ( Better Infection Control )

DYND40708	8Fr (DeLee Tip) Sleeved Catheter, Cup	50/cs
DYND40700	10Fr (Whistle Tip) Sleeved Catheter, Cup	50/cs
DYND40701	12Fr (Whistle Tip) Sleeved Catheter, Cup	50/cs
<b>DYND40702</b>	<b>14Fr (Whistle Tip) Sleeved Catheter, Cup</b>	<b>50/cs</b>



DYND40702

## Open Suction Catheters Kits (Wet)

- Full Range of Sizes ( 6-18Fr)
  - Catheter Kits w/ Solutions
  - Sterile Water ( Preferred )

Item #	Description	Pkg
DYND41476	6Fr, Coiled (DeLee Tip), 2 Gloves, 110ml water	36/cs
DYND41478	Same as above except w/8Fr	36/cs
DYND41470	10Fr, Coiled (Whistle Tip), 2 Gloves, 110ml water	36/cs
DYND41471	Same as above except w/12Fr	36/cs
<b>DYND41472</b>	<b>Same as above except w/14Fr</b>	<b>36/cs</b>
DYND41474	Same as above except w/18Fr	36/cs

### Sterile Saline

DYND41442	14Fr, Coiled (Whistle Tip), 2 Gloves, 110ml saline	36/cs
DYND40442	Same as above with drape	36/cs




DYND41472

**Aloe gloves to help soothe dry hands**



## Open Suction Catheters

Item #	Description	Pkg
<b>- Straight Pack</b>		
DYND41906	6Fr (DeLee Tip)	100/cs
DYND41908	8Fr (DeLee Tip)	100/cs
DYND40160	10Fr (DeLee Tip) No Valve, Calibrated	100/cs
DYND41900	10Fr (Whistle Tip)	100/cs
DYND40990	10Fr (Whistle Tip) Poly-Cath	100/cs
DYND41901	12Fr (Whistle Tip)	100/cs
DYND41162	14Fr (Whistle Tip) No Valve, Calibrated	100/cs
<b>DYND41902</b>	<b>14Fr (Whistle Tip)</b>	<b>100/cs</b>
<b>DYND40992</b>	<b>14Fr (DeLee Tip) Poly-Cath</b>	<b>100/cs</b>
DYND41903	16Fr (Whistle Tip)	100/cs
DYND41904	18Fr (Whistle Tip)	100/cs
DYND40164	18Fr (Whistle Tip) No Valve, Calibrated	100/cs



### - Catheter / Cup & Gloves

DYND40976	6Fr (DeLee Tip), 2 Gloves	50/cs
DYND40978	8Fr (DeLee Tip), 2 Gloves	50/cs
DYND40970	10Fr (Whistle Tip), 2 Gloves	50/cs
DYND40971	12Fr (Whistle Tip), 2 Gloves	50/cs
<b>DYND40303</b>	<b>14Fr (Whistle Tip), 2 Gloves, Mini-Stack</b>	<b>50/cs</b>
<b>DYND40972</b>	<b>14Fr (Whistle Tip), 2 Gloves</b>	<b>50/cs</b>
DYND40973	16Fr (Whistle Tip), 2 Gloves	50/cs
DYND40974	18Fr (Whistle Tip), 2 Gloves	50/cs



DYND40972

## Suction Catheter Tip Descriptions



DeLee Tip

### DeLee Tip

Flat, soft & smooth tip that maximizes newborn, infant and child airways safety. Standard for sizes 6Fr, 8Fr and some 10Fr



Whistle Tip

### Whistle Tip

Angled, soft and smooth tip that minimizes trauma potential to adult large airways and maximizes suctioning capability. Standard in 10Fr, 12Fr, 14Fr, 16Fr and 18Fr sizes.



## Steri-Cath® Closed Suction Systems

### Maxi-Flo® Suction Catheter

Steri-Cath® Systems use the Maxi-Flo® catheter, which is calibrated to assist in proper placement. The soft catheter combines the beveled tip to reduce tracheal and carinal trauma. Each Maxi-Flo® catheter has a large inner diameter for maximum suction efficiency.

### Flexibility without Compromise

Steri-Cath® unibody systems reduce the risk of contamination by eliminating clinician contact with the catheter. It also maintains a closed airway to reduce aerosolization of secretions and oxygen desaturation while allowing constant patient ventilation. These systems are cost-effective and convenient to use.

### Dual Swivel Elbow

An angled patient elbow combined with a dual axis swivel provides torque relief at both the circuit and airway connection (15 mm OD swivel connector at the circuit and 15 mm ID swivel connector at the airway).

### Wiper Seal

An efficient wiper seal maintains PEEP and cleans secretions from the catheter.



PTX610114

PTX615114

### Soft, Flexible Sleeve

A sleeve encloses the catheter to reduce contamination to both health care personnel and patient.

### Single-Step Thumb Valve

The thumb valve allows simple on/off manipulation of the suction source. Raised edges prevent the valve from being inadvertently activated.

Item #	Description	Fr Size	Pkg
<b>Steri-Cath® Single Lumen Suction System</b>			
PTX6100xx	Single lumen, T-piece, swivel & MDI adapter, ET size	12, 14	20/cs
PTX6101xx	Single lumen, T-piece, swivel & MDI adapter, Trach size	12, 14	20/cs
PTX6108xx	Single lumen, T-piece, coude, swivel & MDI adapter, ET size	14	20/cs
PTX6106xx	Single lumen, T-piece, swivel & MDI adapter, 10 each - 15ml saline, ET size	14	20/cs
PTX6118xx	Single lumen, T-piece, swivel & MDI adapter, 3" CORR tubing, ET size	14	20/cs
PTX6119xx	Single lumen, T-piece, swivel & MDI adapter, 3" CORR tubing, Trach size	14	20/cs
PTX6120xx	Single lumen, T-piece, 15x22mm adapter, 3" CORR tubing, Trach size	14	20/cs
<b>Dual Swivel Elbow</b>			
PTX6150xx	Single lumen, dual swivel elbow, one-way valve, ET size	14	20/cs
PTX6151xx	Single lumen, dual swivel elbow, one-way valve, Trach size	12, 14, 16	20/cs

**When ordering, replace "xx" with French size.**

**All Steri-Cath® suction systems include a Trac-Wedge™ and patient label.**

# “Vac-Assist” Suction Aspirator

## HCS7000

Medline’s Vac-Assist Suction Aspirator is simple to use, and has rugged construction to withstand the rigors of homecare use. It is also great for the bedside and crash cart. Features convenient electrical storage. Weighs 12.1 lbs.



-21 inches of mmHg  
Flow>40 liters/minute

AGIBF4053  
Hydrophobic filter  
is standard

Solid base is easy to keep clean  
and features suction cup to keep  
the unit stationary while in use



MIN605  
(45 min of battery  
use at full charge)

Item #	Description	Pkg
HCS7000	Vac-Assist Suction Aspirator (see above)	each
<b>AGIBF4053</b>	Hydrophobic Filter, barbed both ends	3/pk
AGIAG405SET	Tubing and Hydrophobic Filter Set w/Elbow	each
AGIAG405SET3	Tubing and Hydrophobic Filter Set w/Elbow	3/ea
DYND50216	3/16", 6' Sterile Suction Tubing for suction pumps	50/cs
DYND50211	3/16", 20' Sterile Suction Tubing	50/cs
<b>HCSAVF08501</b>	850cc Suction Canister (Float Lid)	12/cs
DYND44600A	800cc Suction Canister (Float Lid) + 20" and 6' Tubing	12/cs
DYND44600	800cc Suction Canister (Filter Lid) + 20" and 6' Tubing	12/cs
<b>MIN605</b>	800 Vacumax Portable Aspirator, AC/DC Power (45 min Battery Suction)	each
HCS701311	Regulator Assembly	each
AGIDF16510	Exhaust muffler filter set (HCS7000)	10/pk



## Gomco Stand Mounted Aspirators and Pumps Dependability and Mobility in One

Model 6037 pump provides suction for all routine suction procedures with one versatile mobile vacuum source.

The ease of operation is simple and consists of four basic steps

1. Set the Mode Switch – Constant/Intermittent
2. Set the Intermittent timing – if needed
3. Turn on the power
4. Adjust the vacuum level

Performance Specifications:

Mode	Vacuum Range	On/Off Cycle Time
Constant	0-250mmHg	On-5-90 seconds
Intermittent	0-200mmHg	Off-23 seconds



GMC6037

Item #	Description	Pkg
GMC6037	Disposable bottle, diaphragm compressor, 1/8 HP cont duty motor, 33"Hx17"Dx23"W, 78 lbs.	each

Gomco suction systems are equipped with an easy-to-read gauge (inches Hg and mmHg), and a suction regulating valve. These models also include illuminated on/off switches.

Model 3021 pump's diaphragm unit is designed for general suction use. Its powerful 1/8 HP motor creates up to 22 in. Hg vacuum

Item #	Description	Pkg
GMC3021	Disposable bottle, 0 to 22 in. Hg vacuum range, 33"Hx17"Dx18"W, 69 lbs.	each

Model 6000 pump is designed to help assure gentle suction for such specialized uses as gastric lavage, abdominal decompression, duodenal and fistula drainage and drainage following prostatectomy.

Item #	Description	Pkg
GMC6000	Glass bottle, 90mm and 120mm Hg vacuum setting, 33"Hx17"Dx18"W, 58 lbs.	each
GMC6002	Disposable bottle, 90mm and 120mm Hg vacuum setting, 33"Hx17"Dx18"W, 58 lbs.	each
GMC2701	Glass bottle, 90mm and 120mm Hg – COMPACT DRAINAGE	each


Model 3040/4042 pumps offer dependable, easy-to-regulate suction for general hospital use. Units feature an insulated, hidden motor and a durable stainless steel top. Model 3041 features a 1/8 HP diaphragm motor, while Model 4042 utilizes a 1/6 HP rotary compressor motor

Item #	Description	Pkg
GMC3040	Disposable bottle, 0 to 22 in. Hg vacuum range, 33"Hx17"Dx18"W, 70 lbs.	each
GMC4042	Disposable bottle, 0 to 25 in. Hg vacuum range, 33"Hx17"Dx18"W, 76 lbs.	each

## Medline Yankauers are Safe and Effective

- Slip-resistant handle for secure grip.
- Smooth inner surface for maximum, non-occluding flow.
- Clear, transparent construction allows for easy assessment of aspirate.
- Shatterproof construction aids patient safety.
- Ribbed connector for various sizes of connecting tubing.
- Optional control vent offered for suctioning control.

Item #	Description	Pkg
<b>Bulb Tip Yankauers, Sterile</b>		
DYND50130	Without control vent, rigid	50/cs
<b>DYND52130</b>	<b>Without control vent, flexible</b>	<b>50/cs</b>
DYND50132	With control vent, rigid	50/cs
DYND52132	With control vent, flexible	50/cs
<b>Flange Tip Yankauers, Sterile</b>		
<b>DYND50140</b>	<b>Without control vent, rigid</b>	<b>50/cs</b>
DYND50142	With control vent, rigid	50/cs
<b>DYND50143</b>	<b>Standard capacity, no vent, flexible</b>	<b>50/cs</b>
<b>DYND50144</b>	<b>Fine capacity, no vent, flexible</b>	<b>50/cs</b>
<b>Bulb Tip Yankauers with Tubing, Sterile</b>		
DYND50135	6' Tubing, without control vent	20/cs
DYND50137	6' Tubing, with control vent	20/cs
<b>DYND50138</b>	<b>10' Tubing, without control vent</b>	<b>20/cs</b>



## Sterile DeLee Mucus Trap with Contro-Vac™ Valve

Transparent 20 cc mucus trap is graduated in 2 cc increments for easy measurement. Extra snap-on cap allows for safe transportation. DeLee catheter has two side eyes, a frosted finish for smooth insertion, and a forward-placed Contro-Vac valve to minimize contact with thumb and fluids. Single-loop flexible packaging with easy-peel lid minimizes storage space.

Item #	Description	Pkg
DYND44108	8Fr Catheter	50/cs
<b>DYND44110</b>	<b>10Fr Catheter</b>	<b>50/cs</b>



DYND44110

## Mucus Specimen Trap

Sterile plastic vials are graduated in 1cc increments. Includes screw-on cap w/latex-free tubing and ID label. (use with any size catheter)

Item #	Description	Pkg
<b>DYND44140</b>	<b>40cc capacity</b>	<b>50/cs</b>
DYND44180	80cc capacity	50/cs



DYND44140

## Use Poole Tips for Gentle Abdominal Suctioning or Spot Suctioning During Surgery

- Rigid construction aids in retraction
- Beveled sleeve openings afford gentle, abdominal suctioning.
- Snap-fit construction permits easy removal of outer sleeve for spot suctioning.
- Two-part tube locks securely for use as a single unit.

Item #	Description	Pkg
<b>Poole Suctioning Instruments, Sterile</b>		
DYND50180	Poole tip	50/cs
<b>DYND50172</b>	<b>Poole tip with 72" non-conductive vinyl tubing</b>	<b>20/cs</b>



DYND50172

## "Patient to the Wall Suction System"



\*DISS Lid top view




Cross-section of canister lid filter



Ask about our Canister Accessories Consignment Program

## Suction Canisters

The filter design in our disposable collection canisters provides significant performance gains. Canister lid features integrated hooks for tubing securement, pre-connected port caps, and a locking lid for maximum safety. \*Unique DISS (screw-on) connection accessory allows canister to attach directly to most suction canisters.



### Stem Lid\*

(For large distances between regulator & canister lid)

Item #	Description	Pkg
GMC20080003	1500cc	48/cs
GMC20080002	2400cc	36/cs
DYND44600	850cc +6' & 20" Tubing	50/cs
DYND44600A	850cc +6' & 20" Tubing, Float Lid	12/cs
DYND44703	1500cc +6' & 20" Tubing	48/cs
DYND44705	1500cc +6' & 20" Tubing +1500cc Safesorb	16/cs

### DISS (Screw-Connect) Lid\*

(Direct to Regulator)

Item #	Description	Pkg
GMC20080004	1500cc	48/cs
GMC20080001	2400cc	36/cs
DYND44805	2400cc +6' & 20" Tubing +3000cc Safesorb	12/cs

\*All have filtered lid except where indicated.



## Medline/H2OR Absorbent Products are the industry's most complete line of premier absorbent products designed to keep the O.R. floor dry, safe and clean.

Create a safe environment in the O.R. by using the Medline/H2OR line of absorbent products to keep the floor free of irrigating fluids. Medline JumboJug, pads and wicks are cost-effective alternatives and competitive products. From floor aspirators, absorbent wicks and pads, to the final fluid collection, Medline/H2OR products absorb fluids on contact to prevent cross-contamination.

### JumboJug Large Volume Suction Canister

- **Sleek profile/integrated handle** reduces the risk of fluid contact and makes transporting simple and efficient.
- **Secure cap** is a snap to use and will accept three devices, including the PowerWick XL.
- **Jumbo collection capacity** 13,000ml (13 Liter) capacity.
- **Fits a standard kick bucket stand**
- **Saves surgical time** by eliminating frequent canister changes when fluid is abundant.

Item #	Description	Pkg
DYNDJJ13K	JumboJug, 13,000 ml	6/bx



Multi-Connection



JumboJug - DYNDJJ13K

### PowerWickXL does not simply absorb – it inhales

This floor aspirator was especially designed for procedures requiring the absorption of a large amount of pooled fluid.

- **Available in 5 foot length.**
- **Low-profile design** is quiet, yet has tremendous vacuum force and pulls in unlimited amounts of fluid in seconds from all sides.
- **Connect to any OR wall vacuum** via JumboJug or standard suction canister.

Item #	Description	Pkg
DYNDPWXL15	PowerWickXL, 5 ft	15/bx



Kick Bucket

PowerWickXL - DYNDPWXL15

# Introducing...

## Medline Lipo-Suction Clear Vinyl Tubing



DYND50110

"We trialed Medline lipo-suction tubing in three different cases today with three different plastic surgeons. . . The physicians and the nurses liked it better than the "special" order tubing we are using because it is longer, more flexible and easier to use. We will now use Medline's lipo-suction tubing."

HealthSouth Outpatient Surgery Center  
Birmingham, Alabama

We manufacture our heavy-duty lipo-suction tubing with 1/4" wall thickness for strength and flexibility.

Item #	Description	Pkg
DYND50110	3/8" x 5/8" x 10'	10/cs
DYND50112	3/8" x 5/8" x 12'	10/cs

- Universal Fit
- Highly Durable
- Easy-to-Use Flexibility
- Utilize with pressures up to 29.5 in Hg (Mercury)

## Sterile Non-Conductive Connecting Tubing

- Available in 5mm (3/16") or 6mm (1/4") inner diameter tubing.
- Ends have female connectors to securely fit any male connector.
- Utilize with pressures up to 22 in Hg (Mercury).

Item #	Description	Pkg
<b>DYND50211</b>	<b>3/16" x 20"</b>	<b>50/cs</b>
DYND50216	3/16" x 6', Male Connector Included	50/cs
<b>DYND50221</b>	<b>3/16" x 10', Male Connector Included</b>	<b>50/cs</b>
DYND50223	3/16" x 12', Male Connector Included	20/cs
DYND50226	3/16" x 20'	20/cs
DYND50240	3/16" x 20" & 6', Tandem Tube Kit, Male Connector Included	50/cs
DYND50246	1/4" x 6', Male Connector Included	50/cs
DYND50251	1/4" x 10', Male Connector Included	50/cs
DYND50252	1/4" x 12', Male Connector Included	20/cs
DYND50253	1/4" x 20', Male Connector Included	20/cs
DYND50230	1/4" x 100', Non Sterile, Bubble Connector Every 3' Bulk (No Female Connector)	1/bx



DYND50221



DYND50211

## Introducing...

### Suction Regulator's Color-Coded, Glow-in-the Dark Gauge Face Provides Easy Readability Under Any Conditions.

- New Soft Touch knobs
- Can be sterilized internally with cold disinfectant
- Weight bearing brass post prevents case from cracking.
- Various inlet and outlet fittings ensure compatibility with existing fixtures in any medical facility
- Five-year warranty
- Minimum maintenance
- Impact-resistant internal components
- Wide variety of inlet and outlet fittings

Item #	Description	Pkg
<b>Continuous/Intermittent Units (0-300mmHg):</b>		
GT881VR30000A	Continuous/Intermittent Regulator, 1 1/4"x1/8NPT inlet, DISS Male outlet	ea
GT881VR300OHC	Continuous/Intermittent Regulator, Ohmeda QC inlet, DISS outlet	ea
GT881VR300OHD	Continuous/Intermittent Regulator, Ohmeda QC inlet, Barb Stem outlet	ea
GT881VR300CHC	Continuous/Intermittent Regulator, Chemetron QC inlet, DISS Male outlet	ea
GT881VR300DHC	Continuous/Intermittent Regulator, Diss FM inlet, DISS Male outlet	ea
<b>Continuous Units (0-300mmHg):</b>		
GT882VR300OHC	Continuous Regulator, Ohmeda QC inlet, DISS outlet	ea
GT882VR300CHC	Continuous Regulator, Chemetron QC inlet, DISS Male outlet	ea
<b>Surgical Units (0-760mmHg):</b>		
GT882VR760OHC	Surgical Regulator, Ohmeda QC inlet, DISS outlet	ea
GT882VR760CHC	Surgical Regulator, Chemetron QC inlet, DISS Male outlet	ea
GT882VR760OHD	Surgical Regulator, Ohmeda QC inlet, Barb Stem outlet	ea



0-160 mmHg  
Gauge



0-300 mmHg  
Gauge



0-760 mmHg  
Gauge

### Medline Sterile Solutions Complement Our Suctioning Line

Pre-packaged sterile saline and sterile water in convenient, ready-to-use containers provide just the right amount of solution necessary for the procedure. Reduces the risk of cross-contamination and potential waste of bulk containers. Bottles are great for device irrigation and suctioning, while the cup size provides the proper amount for suctioning.

Item #	Description	Pkg
PCS1550	250 ml Bottle, Sterile Water, Screw top	24/cs
PCS1650	250 ml Bottle, 0.9% Sterile Saline, Screw top	24/cs
DYND40570	110 ml Cup, Sterile Water	48/cs
DYND40540	110 ml Cup, 0.9% Sterile Saline	48/cs





## SafeSorb™ Fluid Solidifier Is Simply The Safest

SafeSorb is the safest fluid solidifier because it is designed to be poured into an empty suction canister, before a procedure is started★.

Pouring SafeSorb into an empty suction canister offers significant benefits to hospital personnel:

**Safety!** It is safer to pour a solidifier into an empty suction before a procedure starts, rather than opening a port on the suction canister lid to add a solidifier at the end of the procedure, when the suction canister is filled with liquid waste. This benefit:

- **Prevents spilling** the canister contents while pouring the solidifier into a suction canister that is filled with bio-hazardous waste.
- **Eliminates employee exposure** to aerosolized bloodborne pathogens by not having to open a port on the canister lid when the canister is filled with bio-hazardous waste; and,
- **Reduces O.R. turnover time.** With SafeSorb in the suction canister at the start of the procedure, the canister contents are solidified when the procedure is done. This eliminates having to find a bottle of solidifier after the procedure, opening the canister lid and pouring the solidifier into the contaminated canister, then waiting for the canister contents to solidify.

★ SafeSorb can also be added to a suction canister at the conclusion of the procedure. SafeSorb will not materially impact the fluid volume indicated on the graduated lines of a suction canister.

Visit [www.safesorb.com](http://www.safesorb.com):

- For more information about SafeSorb.
- Links to OSHA and Department of Transportation guidelines.
- To quickly and accurately compare the cost of using SafeSorb versus Dornoch, Neptune or a solidifier disinfectant.

Item #	Description	Pkg
DYNSB500CC	Solidifies up to 800cc's of fluid	96/cs
DYNSB1500CC	Solidifies up to 2000cc's of fluid	96/cs
DYNSB3000CC	Solidifies up to 3000cc's of fluid	96/cs
DYNSB1500SP	Solidifies up to 2000cc's of fluid - spout lid for use with canisters having narrow lid ports	96/cs
DYNSB16000CC	Solidifies up to 16,000cc's of fluid - spout lid for use with all canisters and spills	16/cs
DYNSB800T	Dissolvable paper tube solidifies up to 800cc's of fluid	96/cs
DYNSB1500T	Dissolvable paper tube solidifies up to 1500cc's of fluid	48/cs



## Isolyser LTS-PLUS® Fluid Disinfectant/Solidifier

Isolyser LTS-PLUS is EPA registered for the treatment of liquid medical waste that allows disposal directly to a landfill as “**white bag**” waste in 36 states. This allows you not to have to segregate your “white” and “red” bag waste and minimizes exposure.

The patented SAF-T-TOP™ system reduces the need to open a suction canister’s pour spout, thereby, preventing powder spills, aerosolizing, and infectious splashes.

LTS-PLUS utilizes a chlorine-based disinfectant that is preferred for its quick action and non-offensive aroma.

It is available in unit dose vials for treatment of canister or liner liquid medical waste.

Item #	Description	Pkg
MRSLTSP1500	Solidifies 1500cc’s of fluid	100/cs
MRSLTSP2000	Solidifies 2000cc’s of fluid	100/cs
MRSLTSP3000	Solidifies 3000cc’s of fluid	100/cs



### SHAKE WELL



### 5 EASY STEPS!

#### 1. PULL

Leave vacuum on. Pull the patient tubing off the suction canister.

#### 2. PULL TAB

Remove red pull-tab from LTS-PLUS bottle

#### 3. PLACE

Place the LTS-PLUS SAF-T-TOP down over the patient/ortho port on top of the suction canister.

#### 4. PRESS

Press down and allow sufficient time for LTS-PLUS to solidify the waste, because solids will not spill, aerosolize or splash.


#### 5. PROCESS

Process in accordance with local and state regulations.

## Carb-O-Lime Absorbent Granules Efficiently Remove Carbon Dioxide

CarboLime is a granular soda-lime absorbent for the efficient removal of carbon dioxide from closed and semi-closed patient breathing circuits. Contains no potassium hydroxide. Minimum dust levels with the benefits of high surface area and graded particle size. Low bulk density. Low risk of carbon monoxide formation due to good resistance of dry gas desiccation. Low odor due to reliable control of indicator dye concentration.

Item #	Description	Pkg
GMC55010005	Pre-packaged cartridge, 2.5 lb (1.13 kg)	12/cs
GMC55010006	Bag, 3 lb. (1.36 kg)	12/cs
GMC55010003	5-gallon pail	each





## Tracheostomy Care and Cleaning

As a market leader, Medline offers a tracheostomy care line that is unequalled in the U.S. health care system today. We offer products for a wide variety of adult and pediatric needs. This includes our economy, basic and full-component kits and also sterile solutions for catheter cleansing and irrigation of wounds.

Tray Type •Most widely used	Basic			Disposable Cannula	With Solutions		With Suction Catheter	Pediatric Daily Care
Contents	DYND 40582	DYND 40622	DYND 40610	DYND 40620	DYND 40589	DYND 40615	DYND 40580	DYND 40650
Tray w/Tyvek lid	3 comp.	1 comp.	3 comp.	1 comp.	3 comp.	3 comp.	3 comp.	3 comp.
Gloves	2		2	2	2	2	2	2
Drape	1		1	1	1	1	1	
Pipe Cleaners		2	2					
Avant® Gauze (non-woven) 4x4	4	2	2	2	4	2	2	2
Cotton Tip Applicators	2	2	5	5	2	5	2	2
Pliable Cleaning Brush	1	1	1		1	1	1	
Twill Tape 34"	1	1	1	1	1	1	1	1
Avant® Gauze (non-woven) 6 ply Trach-dressing	1 4x4	1 4x4	1 4x4	1 4x4	1 4x4	1 4x4	1 4x4	1 2x2
Hydrogen Peroxide, 1 oz.					2	2		1
Saline 110 ml					1			
Water 110 ml								1
Suction Catheter 14 Fr							1	
Packaging (Qty/cs)	20	30	20	20	20	20	20	20

**TOP  
SELLER!**



DYND40582

### Sterile Tracheostomy Clean & Care Trays

Most sterile sets include a triple compartment tray, moisture-proof drape, pair of vinyl gloves, one premium trach dressing, two cotton-tip applicators, premium gauze sponges, twill tape (1/2" x 34"), and pliable cleaning brush. Peel-down Tyvek lid protects sterile contents during storage and handling. Stacks neatly on shelf or cart.

Item #	Description	Pkg
DYND40582	Features Four Gauze Sponges and Two Cotton-Tipped Applicators	20/cs



## Sterile Tracheostomy Clean & Care Tray with Peroxide and Saline

Complete tracheostomy clean and care tray includes two 1 oz. peel-open foil packs of hydrogen peroxide and one 110 ml cup of saline for convenience and infection control. Triple compartment tray also contains a pair of vinyl gloves; cleaning brush; moisture-proof drape; three 4" x 4", gauze sponges; two cotton-tip applicators; twill tape (1/2"x34"); and a soft, non-linting, non-adherent tracheostomy dressing.

Item #	Description	Pkg
DYND40589	Tracheostomy Care Tray with Peroxide	20/cs



**All-in-one Kit!**



## Medicare Reimbursible A4625 or A4629

### Mini Tracheostomy Clean & Care Tray (Non-Sterile)

Contains the necessary components to perform minor (established) trach cleaning. Includes two 4" x 4", premium gauze sponges; two cotton-tip applicators; two pipe cleaners; trach cleaning brush; and twill tape (1/2"x34") in a non-sterile mini-tray.

Item #	Description	Pkg
DYND40622	Mini Tracheostomy Care Tray	30/cs

## Pediatric-Daily Care Tracheostomy Clean & Care Tray

Ideal for daily wound care of patients requiring a 2" x 2" tracheostomy dressing. Complete tracheostomy clean and care tray includes one 1 oz. peel-open foil pack of hydrogen peroxide and one 110 ml cup of water for convenience and infection control. Triple compartment tray also contains a pair of premium vinyl gloves; two 4" x 4", gauze sponges; two cotton-tip applicators; and a soft 2" x 2" non-linting, non-adherent tracheostomy dressing.

Item #	Description	Pkg
DYND40650	Pediatric Daily Care Tracheostomy Care Tray	20/cs



**Medline Exclusive Tray**

All trays are **latexfree** 



## TrachVox™ Innovative Technology for Tracheostomy Patients




TrachVox is an innovative solution for providing spontaneously breathing tracheostomy patients with the ability to verbally communicate. TrachVox is also a tracheal-connected HME (Heat and Moisture Exchanger). The HME can be depressed which changes the airway flow and allows the patient the ability to speak even with a tracheostomy device.

- Heats and humidifies inspired air of spontaneously breathing patients
- Offers supplemental oxygen through integrated port
- Allows for forced inspiration, expiration and suctioning
- Tested to ISO Standard 9360-1 specifications
- Standard 15 mm connector
- Compact and lightweight
- Patient comfort and minimal trach tube stress

Item #	Description	Pkg
EMS3131010	TrachVox	30

Oxygen connection requires the standard tubing connection to be cut off and the tubing connected to the oxygen inlet port.

**Helps tracheostomy patients become active and maintain a healthy lifestyle.**



**Medicare Reimbursable A7507**

## Tracheostomy Tube Holder (Infant to Adult)

Keeps tracheostomy tube securely in place. Soft strap features an adjustable Velcro closure. Soft compliant foam laminate with elastic adjustment increases patient comfort and accommodates coughing and edema while still providing a snug fit. Soft twill ribbon between collar and tube helps prevent sliding, patient irritation, and helps prevent neck band from becoming soiled.

Item #	Description	Pkg
<b>ZPP501</b>	Tracheostomy tube holder, universal size, 2 piece	20/bx
<b>ZPP301</b>	Tracheostomy tube holder, universal size, 1 piece	20/bx
<b>ZPP501P</b>	Tracheostomy tube holder, pediatric, 2 piece	20/bx
<b>ZPP301P</b>	Tracheostomy tube holder, pediatric, 1 piece	20/bx

Designed as an aid in preventing accidental ventilator disconnections at the trach site and stabilize the trach tube

Item #	Description	Pkg
<b>ZPP401</b>	Vent tie, Anti-disconnect	20/bx

## Endo Tie Endotrach Tube Securing Device

Endo Tie secures using a full length cloth adhesive tape. The cushioned foam pad helps relieve pressure and will not absorb water or adhere to patient's hair.

Item #	Description	Pkg
<b>ZPP707</b>	Endo Tie Trach Tube Securing Device	100/cs



# The *Simplicity* CPAP Line

## HCS903

**The CPAP Product Line demonstrates our commitment to advancements in CPAP technology.**

RestEZ II is engineered with an internal pressure transducer that provides real-time measurements to ensure that the prescribed pressure settings are always being delivered to the patient. The transducer will adjust the blower when a loss of pressure occurs due to resistance in the interface circuit or changes in altitude. This Dynamic Calibration feature can also be used to overcome the resistance inherent in nasal prong interfaces. Built in Auto Compliance feature monitors usage and provides an early identification which creates an opportunity for intervention in the event of a non-compliant patient. The RestEZ Compliance model (HCS903) has enhanced features for measuring and downloading patient compliance data and an Interface Alert, which notifies the patient every 90 days that their interfaces needs replacing.



HCS903 shown with heated humidifier (HCS910)

### HCS903 Specifications

Pressure Range: 3.5 to 20 cm H<sub>2</sub>O  
 Weight, Flow Generator Only: <5.0 lbs (2.27 kg)  
 Warranty: 3 Year warranty  
 Size: 4.7"H x 8.6"W x 9.4"D (11.9cm x 21.8cm x 23.9cm)  
 Power: Input range 100-240V, 50-60 Hz, <120VA  
 (max. power consumption) or 12VDC

### HCS910 Heated Humidifier (Optional)

- Uses only one power supply
- Adjustable temperature dial
- Dishwasher safe water chamber included
- Recommended to boost compliance

## BiPAP Plus M-Series



Item #	Description	Pkg
HDTDS600	BiPAP Plus	each
HDT102257	BiPAP Heated Humidifier	each





SEE50725

## CPAP Accessories

Item #	Description	Pkg
<b>SEE50725</b>	<b>I.Q. adjustable gel mask</b>	<b>each</b>
<b>AEIHRML</b>	Headrest with Nasal Pillows	each
AGIAC302328	Simple strap headgear for SEE50725	each
AGICP11602SL	CPAP mask, medium with vented swivel	each
AGICP11602SA	CPAP mask, small with vented swivel	each
AGIHCG140C	18" CPAP tubing	each
AGICP11524	Vented swivel elbow adapter for all contour masks	each
<b>HCSCM003X</b>	Mojo Full Face Mask, Incl. - Headgear -(X= S,M, L)	each
	Gel - A7034/7035	
<b>HCSCM007X</b>	CPAP/BiPaP FF Mask, Incl. - Headgear -(X= S,M, L)	each
	Silicone -A7034/7035	
<b>HCSCM008X</b>	CPAP/BiPaP Nasal Msk, Incl. - Headgear -(X=S,M, L)	each
	Silicone -A7034/7035	



HCSCM003X - Mojo



AEIHRML - Nasal Prongs



HCSCM007X - Reusable



HCS2367133B (Display)

## CPAP Accessories Cleaner

**Purpose:** Maximize Patient CPAP Therapy Compliance

- Alcohol-Free and Solvent-Free ( Natural )
- No odor or residue left behind after cleaning
- Cleaner dries on its own and does not need to be rinsed ( no rinsing )
- Wipes will last 2 months for every patient, if they use one daily
- The 1.5oz & 8oz Traveler Spray canisters are propelled by Compressed Gas

### Benefits:

Allows the mask usage to be maximized ( up to 3 months )

Minimize extra mask expense without reimbursement ( 3 months – CMS )

Allows the mask to fit properly without over-tightening


Maximize profitability with fewer service calls and expense



HCS635871164



HCS635871165



HCS635871639

HPCS#	Item #	Description	Pkg
<b>No code</b>	HCS2367133B	Medline Store Display (restock what gets sold) Includes: 9 ea CPAP Wipes, 4 ea 32 Oz Concentrate 3 ea 8 Oz Spray 2ea 1.5 Oz Spray	each
<b>No code</b>	HCS635871639	CPAP Wipes	12/cs
<b>No code</b>	HCS635871164	1.5 Oz Spray	24/cs
<b>No code</b>	HCS635871165	8 Oz Spray	12/cs
<b>No code</b>	HCS635871849	32 Oz Concentrate (tubing/headstrap)	4/ea

## CPAP Headgear

### Simple Strap

The Medline Simple Strap Headgear can be used for many of the most common masks on the market. The convenience of the Simple Strap is unsurpassed because one size fits all and its simple design is the ultimate user friendly headgear. The vented head piece minimizes perspiration associated with other types of headgear. It is available in blue or black.

Item #	Description	Pkg
AGIAC302328	Simple Strap, Universal, Blue	each
AGIAC302328B	Simple Strap, Universal, Black	each



AGIAC302328B

### Coolcap

The Medline Coolcap is designed for the more active patient whom has difficulty with other headgear. It is designed to minimize perspiration due to patient heat retention. The Medline Coolcap will fit most any manufacturer's most popular mask for your convenience. Available in small, medium and large sizes.

Item #	Description	Pkg
AGIAC302458	Coolcap, Medium	each
AGIAC302439	Coolcap, Large	each
AGIAC302442	Coolcap, Extra-Large	each



AGIAC302439

For additional options - Contact your Medline Respiratory Product Support



AGIHCG142C

## CPAP Tubing

Our Medical bore tubing is manufactured of the highest quality to meet your demanding needs. Tubing features a smooth interior which allows for excellent air flow, and 22mm rubber cuffs or integral cuffs. Tubing maintains cross sectional diameter at bends.

Item #	Description	Pkg	
<b>HCSSBTCAP</b>	6', cream	10/cs	<b>Best Value!</b>
AGIHCG145C	8', cream	each	
AGIHCG146C	10', cream	each	
<b>AGIHCG142C</b>	4', grey	each	
<b>AGIDS143C</b>	72" Economy CPAP circuit	each	



AGIDS143C - Economy tubing





## Introducing the AutoVent 4000 Series Transport Ventilator with optional CPAP

- Pneumatically powered....durable and light, with all basic models weighing less than 5 lbs.
- Independent settings for the critical ventilation parameters of Breaths Per Minute, Inspiratory Time, Tidal Volume and Airway Pressure Relief.
- Manual Breath Button and Anti-Breath Stacking standard on all models.
- Large color-coded Manometer Gauge along with a full set of alarms, to ensure optimal performance and patient safety.
- Optional Air Mixing provides 65% O<sub>2</sub> concentration for extended cylinder life during lengthy transports.
- The capabilities of ventilation and CPAP (5-20 cwp) in one system
- Electronic Alarm package available





## L760

- Pneumatic Operation
- Adjustable Tidal Volume
- Adjustable BPM
- Child and Adult Settings
- Internal Demand Valve
- Anti-Breath Stacking
- Manometer Readout
- Adjustable Pressure Relief with Alarm
- Low Source Gas Pressure Alarm
- Color-coded Labels and Gauge for initial Settings
- Manual Breath Button

## L761

- Pneumatic Operation
- Adjustable Tidal Volume
- Adjustable BPM
- Child and Adult Settings
- Internal Demand Valve
- Anti-Breath Stacking
- Manometer Readout
- Adjustable Pressure Relief with Alarm
- Low Source Gas Pressure Alarm
- Color-coded Labels and Gauge for initial Settings
- Manual Breath Button
- **Air Mix Mode**

## L762

- Pneumatic Operation
- Adjustable Tidal Volume
- Adjustable BPM
- Child and Adult Settings
- Internal Demand Valve
- Anti-Breath Stacking
- Manometer Readout
- Adjustable Pressure Relief with Alarm
- Low Source Gas Pressure Alarm
- Color-coded Labels and Gauge for initial Settings
- Manual Breath Button
- **Air Mix Mode**
- **Electronic Alarm Package**

## L763

- Pneumatic Operation
- Adjustable Tidal Volume
- Adjustable BPM
- Child and Adult Settings
- Internal Demand Valve
- Anti-Breath Stacking
- Manometer Readout
- Adjustable Pressure Relief with Alarm
- Low Source Gas Pressure Alarm
- Color-coded Labels and Gauge for initial Settings
- **Electronic Alarm Package**

## All models available with optional CPAP

Item #	Description	Pkg
<b>AutoVent</b>		
L760	AutoVent 4000	each
L760CPAP	AutoVent 4000 with CPAP	each
L761	AutoVent 4000 with Air Mix Mode	each
L761CPAP	AutoVent 4000 with Air Mix Mode and CPAP	each
L762	Autovent 4000 with Air Mix Mode and Electronic Alarm Package	each
L762CPAP	Autovent 4000 with Air Mix Mode and Electronic Alarm Package with CPAP	each
L763	AutoVent 4000 with Electronic Alarm Package	each
L763CPAP	AutoVent 4000 with Electronic Alarm Package and CPAP	each
<b>Accessories</b>		
GMCL599130	Ventilator Circuit, 3', Swivel, Adult #5 Mask	10/cs
GMCL535026	Oxygen Hose, 6' with DISS fittings	each
GMCL270220	Oxygen Regulator, 2 DISS, 1 Barb outlet, 870CGA	each
GMC800149	Internal Air Filter (replace when soiled)	3/pkg
GMC800148	CPAP Ventilator Circuit 3', EXH Valve, Cuffed, Adult #5 Mask	6/cs

**All systems comply with International Standards of FDA and Canada**

### Specifications:

#### L760 AutoVent 4000

Dimensions: 2.9" x 9.6" x 6.5"

Weight: 3.7 lbs.

Weight w/CPAP: 5.2 lbs.

#### L761 AutoVent 4000 with Air Mix Mode

Dimensions: 2.9" x 9.6" x 6.5"

Weight: 4.2 lbs.

Weight w/CPAP: 5.7 lbs.

#### L762 Autovent 4000 with Air Mix Mode & Electronic Alarm Package

Dimensions: 2.9" x 9.6" x 8.0"

Weight: 4.9 lbs.

Weight w/CPAP: 6.4 lbs.

#### L763 AutoVent 4000 with Electronic Alarm Package

Dimensions: 2.9" x 9.6" x 8.0"

Weight: 4.4 lbs.

Weight w/CPAP: 5.9 lbs.



## Introducing the EPV100 Portable Ventilator for Mass Casualty Preparedness

- Designed to Meet Surge Capacity
- Light and Easy to Operate
- Powered by Compressed Gas
- Suitable for Stockpiling or Everyday Use.



Item #	Description	Pkg
EPV100	Portable Ventilator, 6' disposable O2 tubing, 3' ventilator circuit with swivel, adult disposable cuffed mask, filter	each
<b>Accessories</b>		
L599-130	3' Ventilator Circuit, with swivel, adult disposable cuffed mask, filter	10/cs
L535026	6' Oxygen Hose, reusable	each
L270-220	Oxygen Pressure Regulator, 0-25 lpm, 2 check valves, hose barb outlet*	each

\*Additional LSP Oxygen Pressure Regulators with check valves are available

## iVent201 Making Patient Care Simpler

Use it in sub-acute and nursing homes.

Easy software upgrades.

Make your patients more comfortable than ever.

Easy to Learn Easy to Use

The iVent201's Adaptive Flow™ and Easy Exhale™ features provide unprecedented patient comfort. These breakthrough technologies keep patients from "fighting the ventilator" by letting their needs spontaneously determine flow rate.



Item #	Description	Pkg
VRS501I2201HC	<b>iVent HC</b> (Home Care/Nursing Home): Modes:Volume Control (A/C or SIMV), CPAP/PSV Low Pressure Oxygen.	each
VRS501I1201IC	<b>iVent201 IC</b> (Intensive Care): Modes: Pressure Control (A/C or SIMV),Volume Control (A/C or SIMV), CPAP/ PSV Software Packages Included: Pressure, Flow and Volume Waveforms,Trending and Respiratory Diagnostics, Internal Oxygen Mixer with Sensor, High and Low Pressure Oxygen, MRI Compatible	each

## HME Filter Performance Data and Specifications

Product Code	Physical Appearance	Connectors (ISO)	Dead Space (cc)	Weight (gm)	Type Of Filtration	BFE %	VFE %	CO <sub>2</sub> Port	*Moisture Output (Mg H <sub>2</sub> O/L)
HCSMC700	Angled	15 mm I.D./22 mm O.D. x 15 mm O.D.	22.5	18.8	Electro-Static	99.99	99.99	Yes	30.4
HCSMC700NF	Angled	15 mm I.D./22 mm O.D. x 15 mm O.D.	22.5	18.8	—	—	—	Yes	30.4
MC-701	Angled	15 mm I.D./22 mm O.D. x 15 mm O.D.	22.5	18.8	Electro-Static	99.99	99.99	No	30.4
HCSMC701NF	Angled	15 mm I.D./22 mm O.D. x 15 mm O.D.	22.5	18.8	—	—	—	No	30.4
HCSMC750	Straight	15 mm I.D./22 mm O.D. x 15 mm O.D.	22.5	18.7	Electro-Static	99.99	99.99	Yes	30.4
MC-750NF	Straight	15 mm I.D./22 mm O.D. x 15 mm O.D.	22.5	18.7	—	—	—	Yes	30.4
HCSMC751	Straight	15 mm I.D./22 mm O.D. x 15 mm O.D.	22.5	18.7	Electro-Static	99.99	99.99	No	30.4
MC-751NF	Straight	15 mm I.D./22 mm O.D. x 15 mm O.D.	22.5	18.7	—	—	—	No	30.4
HCSMC775A	Straight	15 mm I.D./22 mm O.D. x 22 mm I.D.	22.5	18.7	Electro-Static	99.99	99.99	No	30.4
HCSMC775NF	Straight	15 mm I.D./22 mm O.D. x 22 mm I.D.	22.5	18.7	—	—	—	No	30.4
HSK763	Straight	15 mm I.D./22 mm O.D. x 22 mm I.D.	48.5	35.0	Electro-Static	99.99	99.99	No	33.2
HSK761	Straight	15 mm I.D./22 mm O.D. x 22 mm I.D.	48.5	35.0	—	—	—	No	33.2
MC-790	Straight	15 mm I.D. x 15 mm O.D.	2.0	3.7	—	—	—	No	32.1

### Medicomp MC-750 (Adult & Pediatric HME/BFE/VFE) Filter

The sample was tested at a test condition of 0.5 L tidal volume, a frequency of 20 breaths per minute and an I:E ratio of 1:2. The moisture output of the MediComp MC-750 was 30.4 mg H<sub>2</sub>O/L.

## HME Filters

Item #	Description	Pkg
<b>HCSMC700</b>	HME w/bacterial viral filter w/CO <sub>2</sub>	50/cs
HCSMC750	HME w/bacterial viral filter, angled w/CO <sub>2</sub> (22mmID x 15mmOD)	50/cs
<b>HCSMC701NF</b>	HME only (angled)	50/cs
<b>HCSMC751</b>	HME w/filter, No CO <sub>2</sub> Port (15mmOD x 15mmID)	50/cs
HCSMC775A	HME w/filter, No CO <sub>2</sub> Port (15mmOD x 22mmID)	50/cs
HCSMC751NF	HME/ Pressure (15mmOD x 15mmOD)	50/cs
<b>HCSAVF1675</b>	HME Flex/F15/F22	50/cs



HCSAVF1675



## A

Adapters	9
Aeromist Epic Nebulizer/Compressor System	24
Aerosol Drainage Bags	19
Aerosol Masks	19
Aerosol Therapy	17-19
Aerosol Tubing	18
Air Compressor (50 psi)	18
AMBU Bags	6
Asthma Management	20-25

## B

Bubble Humidification	9
BiPAP	43

## C

Cannula	7
Carts	15
Connectors	9
Connecting Tubing	36
Corrugated Aerosol Tubing	18
CPAP Therapy	43-45
CPAP Masks	44
CPAP Tubing	45
CPAP Headgear	44
CPR Mask	5

## D

Digit Finger Oximeter	10
-----------------------	----

## E

EndurO <sub>2</sub> Oxygen Concentrator	4
---	---

## F

Face Shields	5
Filters	49
Fluid Disinfectant	39

## G

Gomco Pumps	32
-------------	----

## H

Handheld Pulse Oximeter	10
Headgear	45
Heat Moisture Exchangers (HME's)	49
Holding Chambers	21

## I

Incentive Spirometer	21
Isolyer LTS Fluid Disinfectant/Solidifier	39

## L

Large Volume Nebulizers (Sterile H <sub>2</sub> O)	18
Leverhaler	21

## M

Manual Resuscitators	5
Masks CPAP	44
CPR	5
Oxygen	8
Trach	19
MaxiCompressor	18
Metered Dose Inhaler Holding Chamber (Leverhaler)	21
Mouth to Mask Resuscitator	5
Mucus Traps	33

## N

Nasal Cannula	7
Nebulizer Compressor	22-24
Nebulizer Kits	25

## O

Open Suction Catheters and Kits	28-29
Opti-Chamber	21
Oxygen Analyzer	5

Oxygen Cannulas	7
Oxygen Carts	15
Oxygen Concentrators	5
Oxygen Connectors and Adapters	9
Oxygen Conservers	12-13
Oxygen Flow Meter/Adapters	16
Oxygen Masks	8
Oxygen Regulators	14
Oxygen Therapy	3-16
Oxygen Tubing	7

## P

Peak Flow Meters	21
Poole Tips	33
Pulse Oximetry	10-11

## R

Regulation Devices	14
--------------------	----

## S

Safesorb Fluid Solidifier	38
Saline Vials	25
Soft-Touch Oxygen Cannula	7
Solidifiers	38-39
Sterile Saline/Water	37
Suction	26-39
Suction Canisters	34
Suction Catheters	27-30
Suction Irrigation Solutions	33
Suction Pumps	31-32
Suction Regulators	37
Suction Tubing	36

## T

Tracheostomy Care and Cleaning	40-42
Tracheostomy Care Trays	40-41
Trach Masks	19
TrachVox	42

## V

Vacumax Suction Aspirator	31
Ventilators	46-48
Ventilation Management	46-49
VixOne Nebulizer Circuits	25

## Y

Yankauers	33
-----------	----

## Quality Respiratory Care Products You Can Always Depend On.

Quality is key in any business and especially where people's lives can depend on the performance of our products, quality is everything. That is why everything we do starts and ends with our customer. Our livelihood rests solely on our ability to provide quality respiratory and anesthesia care products where and when our customers need them. We stand behind each and every one of our products and guarantee they will perform flawlessly and without complication.

One of the strongest demonstrations of Medline's commitment to the production of consistently high quality products is the adherence to the Food and Drug Administration's (FDA) Quality System Regulation (QSR), International Organization for Standardization (ISO) and European Norm (EN) standards. Medline has divisions registered to ISO9001 and EN46001 Quality Management Systems. Bestowed by AFAQ, the French association for third party certification of quality systems, ISO certification is granted only after rigorous analysis of a company's quality systems.

*Process control from equipment  
design through validation enables  
us to provide customers the  
consistency they expect.*



## Medline is Your One-Stop Shop

Medline offers everything you need to become the only place you purchase your respiratory products from. With quality products, excellent customer service, and satisfaction guaranteed; Medline wants to become your one source.



1-800-MEDLINE (633-5463)

[www.medline.com](http://www.medline.com)

©2008 Medline Industries, Inc. Medline is a registered trademark of Medline Industries, Inc. One Medline Place, Mundelein, IL 60060  
LTS-PLUS is a registered trademark and SAF-T-TOP is a trademark of Microtek Medical Holdings, Inc. SafeSorb is a registered trademark of DiSorb Systems, Inc.  
Spiritus is a trademark of Viasys Inc. SleepStrip is registered trademark of SLP Ltd.  
We reserve the right to correct any errors that may occur in this brochure.

Medline Canada 1-800-396-6996

[www.medline.ca](http://www.medline.ca) [canada@medline.com](mailto:canada@medline.com)

Not all products are currently available for sale in Canada.  
Please contact your Account Representative for additional information.

MKT208208/LIT085/ZPN6/10M