

PATIENT SLINGS*

Standard Slings

ITEM #	DESCRIPTION
30137	Standard seat/back sling - Canvas 350lb weight capacity
30138	Standard seat/back sling w/ commode - Canvas 350lb weight capacity
30132	Standard seat/back sling - Nylon 350lb weight capacity
30133	Standard seat/back sling w/ commode - Nylon 350lb weight capacity



*All slings have been independently weight capacity certified by Intertek Testing Services according to ISO10535 standards.

PATIENT SLINGS*

Full Body Mesh Slings

ITEM #	DESCRIPTION
30110 Medium	Mesh full body medium 450lb weight capacity <i>(best fit 100lb-195lb)</i>
30111 Large	Mesh full body large 450lb weight capacity <i>(best fit 150lb-300lb)</i>
30118 XX-Large	Mesh full body XX-large bariatric 600lb weight capacity <i>(best fit 290lb-600lb)</i>



Full Body Mesh Slings W/ Commode

ITEM #	DESCRIPTION
30114 Medium	Mesh full body Medium 450lb w/ commode opening <i>(best fit 100lb-195lb)</i>
30115 Large	Mesh full body Large 450lb w/ commode opening <i>(best fit 150lb-300lb)</i>
30116 X-Large	Mesh full body X-Large 450lb w/ commode opening <i>(best fit 210lb-400lb)</i>
30119 XX-Large	Mesh full body XXL bariatric 600Lb w/ commode opening <i>(best fit 290lb-600lb)</i>



*All slings have been independently weight capacity certified by Intertek Testing Services according to ISO10535 standards.

PATIENT SLINGS*

Full Body Solid Slings

ITEM #	DESCRIPTION
30112 Medium	Solid fabric full body medium 450lb weight capacity (<i>best fit 100lb-195lb</i>)
30113 Large	Solid fabric full body large 450lb weight capacity (<i>best fit 150lb-300lb</i>)
30117 X-Large	Solid fabric full body extra large 450lb weight capacity (<i>best fit 210lb-400lb</i>)



Divided Leg Slings

ITEM #	DESCRIPTION
30100P Petite	Solid divided leg Petite 450lb (<i>best fit 140lb-200lb</i>)
30100 Small	Solid divided leg Small 450lb (<i>best fit 200lb-400lb</i>)
30101 Medium	Solid divided leg Medium 450lb (<i>best fit 200lb-450lb</i>)
30102 Large	Solid divided leg Large 450lb (<i>best fit 200lb-450lb</i>)



*All slings have been independently weight capacity certified by Intertek Testing Services according to ISO10535 standards.

PATIENT SLINGS*

Sit to Stand / Transfer Slings

ITEM #	DESCRIPTION
30130	Sit to Stand Sling - Solid 450lb weight capacity
30131	Sit to Stand Sling – Solid Large 450lb weight capacity
30134	Transfer Sling - Solid 450lb weight capacity
30136	Transfer Sling – Solid Large 450lb weight capacity



*All slings have been independently weight capacity certified by Intertek Testing Services according to ISO10535 standards.