

BLICKMAN



Includes:

- IV Stands
- IV Stand Accessories
- IV Irrigation Safety Systems

IV Stands
Optimal continuous fluid delivery

OUR MISSION

Customer Focus

Exceed customer expectations through a commitment to continuous improvement of our products and services.

Quality Products

Offer the marketplace a broad range of high quality, cost effective, and durable products of superior craftsmanship.

Organizational Excellence

Strive for excellence in all we do by fostering a participative, team based environment, which nurtures the active involvement of all members of the organization.

Marketplace Recognition

Earn the respect, admiration, and loyalty of our people, customers, and suppliers by treating them with dignity and respect, embracing diversity, and developing and maintaining long term partnerships.

Collaborative Workplace

Encourage all employees to work safely and efficiently, contribute innovative ideas, and personally benefit from the ongoing success of the company.

The logo for Blickman, featuring the word "BLICKMAN" in a bold, white, sans-serif font. The text is contained within a dark blue rectangular box with a thin white border. The box is set against a light gray background that has a subtle gradient and a slight shadow effect.


THE BLICKMAN DIFFERENCE

IV Stands

- Our IV poles are optimal for maintaining continuous fluid delivery for every patient, available in a wide range of adjustable heights with a variety of easy locking options that can be operated with one hand.
- Our IV poles are all available in high quality chrome and stainless steel, and set on 4, 5, or 6 leg bases with casters, no projections, and great balance/tip resistance. Some bases are gravity weighted to help compensate for heavier loads. All bases have the ability to spin 360° in tight spaces.
- Most of our models have 2, 4, 6 and 8 (rake) hook adapters available. All come with colored secure grip tips.

Key Product Features

*See page 9 for more details on all optional features.



IV Stands

4 Leg

4 LEG CHROME IV STANDS

Description: We offer two 4 leg chrome IV stand models, one with a standard base, and one with a heavyweight base. Both base types include Wall Saver™ tips to protect against wall scuff marks. They both move on 2" plastic twin wheel hooded swivel casters and they both come standard with two hooks. A four hook removable adapter is available for model 1305.



Height adjustment:

47 1/2" - 85 1/2"

Lock:


Twist

Base:

22" diameter
Apoxy coated
chrome

Casters

2" plastic twin wheel
hooded swivel
casters



Height adjustment

51 1/2" - 93"

Lock:

Tru-loc™

Base:

22" diameter
Chrome plated
heavy bar steel

Casters

2" plastic twin wheel
hooded swivel
casters

Name	Standard Base	Heavyweight Base
# of Hooks	2 Welded Hooks	2 Removable Hooks
Model #	1305	1310
Item #	0561305000	0561310000
Shipping Weight	13 lbs.	14 lbs.

IV Stands

4 Leg

4 LEG STAINLESS STEEL IV STANDS

Description: We offer two types of 4 leg stainless steel IV stands, one with a standard base, and one with a heavyweight base. They come standard with either 2 or 4 removable hooks.
Contact customer service for 6 or 8 (rake) hooks.



Name	Standard Base
# of Hooks	2 Hooks
Model #	1410SS
Item #	0571410000
Shipping Weight	13 lbs.
# of Hooks	4 Hooks
Model #	1410SS-4
Item #	0571410400
Shipping Weight	14 lbs.

Height Adjustment: 51 3/4" - 94"

Mechanism: Tru-Loc™

Base Material: Stainless steel with Wall Savers™, 22" diameter

Casters: 2" hard rubber swivel casters



Name	Heavyweight Base	
Lock Type	Twist Lock	Thumb Control
# of Hooks	2 Hooks	
Model #	8890SS	8889SS
Item #	0528890000	0518889000
Shipping Weight	24 lbs.*	24 lbs.*
# of Hooks	4 Hooks	
Model #	8890SS-4	8889SS-4
Item #	0528890400	0518889400
Shipping Weight	24 lbs.*	24 lbs.*

Height Adjustment: 66" - 100"

Mechanism: Twist Lock, Thumb Control, Foot Operated

Base Material: Heavyweight stainless steel flat bar, 22" diameter.

Casters: 2" roller blade swivel casters

*Shipping weight with an asterisk is the combined shipping weight of separate parts.




IV Stands

5 Leg

5 LEG IV STANDS


Description: Our 5 leg IV stands come in either chrome or stainless steel poles and stand on a 24" diameter, polished die cast aluminum base with Wall Saver™ caps. They come standard with either two or four removable hooks. Contact customer service for 6 or 8 (rake) hooks.



Name		Chrome Pole	
Pole	Heavy Duty (11/4")	Standard (1")	
# of Hooks	2 Hooks		
Caster Heights	3" Rubber	2 1/2" Plastic	
Height Adjustments	56" - 100"	52 1/2" - 93 1/2"	
Model #	1350-2	1315	
Item #	0561350200	0561315000	
Shipping Weight	17 lbs.*	16 lbs.*	
# of Hooks	4 Hooks		
Caster Heights	3" Rubber	2 1/2" Plastic	
Height Adjustments	56" - 100"	52 1/2" - 93 1/2"	
Model #	1350 -4	1315-4	
Item #	0561350400	0561315400	
Shipping Weight	17 lbs.*	16 lbs.*	

Height Adjustment Mechanism: Tru-Loc™
 Base Material: 5 Leg, Polished die cast aluminum, 24" diameter.
 Casters: 3" hard rubber OR 2 1/2" twin wheel swivel casters.





Name		Stainless Steel Pole		
Lock Type	Tru-Loc™	Twist Lock	Thumb Control	
# of Hooks	2 Hooks			
Caster Heights	2" Rubber	3" Rubber	3" Rubber	
Height Adjustments	52" - 94"	74" - 110"	67" - 98"	
Model #	1415SS	7792SS	7794SS	
Item #	0531415000	0537792000	0517794000	
Shipping Weight	16 lbs.*	15 lbs.*	16 lbs.*	
# of Hooks	4 Hooks			
Caster Heights	2" Rubber	3" Rubber	3" Rubber	
Height Adjustments	52" - 94"	74" - 110"	67" - 98"	
Model #	1415SS-4	7792SS-4	7794SS-4	
Item #	0531415400	0537792400	0517794400	
Shipping Weight	16 lbs.*	16 lbs.*	17 lbs.*	

Height Adjustment Mechanism: Tru-Loc™, Twist Lock, Thumb Control
 Base Material: 5 Leg, Polished die cast aluminum, 24" diameter.
 Casters: 2" or 3" hard rubber swivel casters.

*Shipping weight with an asterisk is the combined shipping weight of separate parts.

MR option available, see MR section for more details.



IV Stands

5 Leg

5 LEG WASHABLE IV STAND



Description: Our Hackensack 5 leg IV stand is cart washable as it was built specifically to withstand high temperatures and cleaning fluids with weep holes in the pole and in the base for fluid drainage. This model comes with a stainless steel pole with a twist lock height adjustment mechanism. Choose from 2, 4, and 8 (rake) removable hooks.

Features Include:

- ✓ 22 7/8" diameter, 5 leg black composite base.
- ✓ 2 1/2" twin wheel hooded swivel casters.
- ✓ 1" outside diameter base pole, and 3/4" outside diameter extension pole.
- ✓ Height extends from 74" to 110" off the ground.



Name	Washable IV Stand
Lock Type	Twist Lock
# of Hooks	2 Hooks
Model #	7795SS-2
Item #	0539A02302
Shipping Weight	16 lbs.*
# of Hooks	4 Hooks
Model #	7795SS-4
Item #	0539A02304
Shipping Weight	16 lbs.*
# of Hooks	8 Hooks
Model #	7795SS-RH
Item #	0539A023RH
Shipping Weight	16 lbs.*

IV Stands

6 Leg

6 LEG IV STANDS



Description: Our 6 leg IV stands are designed for hospital applications requiring heavier loads by using low center of gravity bases. They come in either a chrome or a stainless steel pole with a Tru-Loc™ height adjustment mechanism that allows the pole to reach a minimum of 52" high off the ground and a maximum of 90" off the ground. Our 6 leg IV stands come standard with 4, 6, and 8 (rake) removable hooks.

Features Include:

- ✓ 23 3/8" diameter base with low center of gravity for improved stability.
- ✓ Heavy duty 3" roller blade swivel casters for easy maneuverability.
- ✓ Steel powder coated finished base for durability (many colors available).
- ✓ Accommodates up to 4 infusion pumps and up to 8 IV solution bags while preventing tubing line entanglement.



Powder Coated Base

Stainless Steel Base

Name	Chrome Pole	SS Pole	SS Pole
# of Hooks	4 Hooks		
Model #	1370-4	1370SS-4	1371SS-4
Item #	0561370400	0541370400	0541371400
Shipping Weight	29 lbs.*	29 lbs.*	26 lbs.*
# of Hooks	6 Hooks		
Model #	1370-6	1370SS-6	1371SS-6
Item #	0561370600	0541370600	0541371600
Shipping Weight	29 lbs.*	29 lbs.*	26 lbs.*
# of Hooks	8 (rake) Hooks		
Model #	1370-RH	1370SS-RH	1371SS-RH
Item #	0561370RH1	05413700RH	05413710RH
Shipping Weight	29 lbs.*	29 lbs.*	26 lbs.*



Custom Powder Coated Colors Available Upon Request

Please email us at info@blickman.com

IV Stands

Options

IV Hook Configurations



Part #	Description	Shipping Weight
780063700L	2 hooks	2 lbs.
780063800L	4 hooks	2 lbs.
780066900L	6 hooks	2 lbs.
7802808750	8 (rake) hooks	6 lbs.

Lock Options



Thumb Lock

Tru-Lock™

Twist Lock

Base Options

4 Leg Chrome Bases



Standard



Heavyweight

4 Leg Stainless Steel Bases



Standard



Heavyweight

5 Leg Bases



Standard



Washable

6 Leg Bases



Powder Coated



Stainless Steel

IV Stands

Accessories

IV POLE ATTACHMENTS



Wire Basket



Guidance Handle



Power Strip



Accessory Clamp



Drainage Bag Hook



Pump Support Tray



Multi-pump Support



Tank holder

Part #	Description	Shipping Weight
3501389001	Wire Basket	9 lbs.
3501381000	Guidance Handle (MR Option Available)	4 lbs.
9101107600	Power Strip	4 lbs.
7801380000	Accessory Clamp	4 lbs.
7801382000	Drainage Bag Hook	1 lbs.
7801384000	Pump Support Tray	2 lbs.
7801385000	Multi-pump Support	2 lbs.
7801386000	Tank Holder	5 lbs.

MR options available where indicated, see MR section for more details.



IV Stands

Construction Specifications

Construction Specifications

Universal Dimensions	Height Range: 47" – 110" Base Diameter Range: 22" – 24"
Materials	304 stainless steel Chrome
Gauge	Stainless Steel IV Stands: Base Pole: 1" Diameter, 16GA; Extension Pole: ¾" Diameter, 18GA Exception: 1370/1371 series Base Pole: 1 ¼" Diameter, 16 GA; Extension Pole: 1" Diameter, 18GA
Welding	All welding by heliarc method and exposed welds are ground smooth.
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish.
Casters	See page 9.
Notes	MR conditional IV stand option available where indicated.

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Note for: <i>IV Stands</i>	Parts shipped separately to avoid damage and lower costs. Shipping weight with an asterisk is the combined shipping weight of separate parts. Some assembly required.
--------------------------------------	---

THE BLICKMAN DIFFERENCE

IV Irrigation Safety Systems

- Our IV Irrigation Safety Systems (IVTs) offer superior safety, efficiency, and financial value for fluid intensive IV irrigation procedures, such as arthroscopic, hysteroscopic, endoscopic, TURP, etc.
- They maintain continuous fluid delivery by offering infinitely adjustable heights for up to eight 3 liter or 5 liter IV bags on four independent channels.
- Their easy, one handed operation allows nurses to attach IV bags at comfortable heights, starting at 5 ft., and then raise bags up to 8 ft. (or to 9 ft. with custom extender) simply by pushing down on “squeeze/release” handles, an ergonomic solution that meets OSHA standards and greatly reduces the threat of back injuries.

Key Product Features

Up to eight 3 liter or 5 liter IV bags on four independent channels.

High quality, easy maintenance, light weight aluminum.




Locking casters for safe mobility and 360° spin in tight spaces.



Push down “Squeeze/Release” handles to raise IV bags.

Indefinite adjustable heights from 5 ft. to 8 ft. (9 ft. with custom extender).

26 1/2" diameter, inverted, 5 legged base resistant to tipping.



WHY CHOOSE OUR IVT?



Safe:

- Their easy, one handed operation allows nurses to attach IV bags at comfortable heights, simply by pushing down on “squeeze/release” handles, an ergonomic solution that meets OSHA standards and greatly reduces the threat of back injuries.
- They are constructed of hard anodized aluminum with stain resistant finishes. So, no surfaces can chip, flake, or rust to introduce particle matter into a sterile or surgical environment.
- Their lowest positions are cushioned by shock absorbers, thus sparing the common banging known to occur with the standard IV poles, which results in tearing of bag tabs and spilling remaining fluids onto the O.R. floor.
- The steady, five legged bases resist potential tipping hazards, eliminating the need for “sand bagging”.

Efficient:

- There are no electrical aspects that require overnight charging. They are always ready, whenever you are ready.
- You can stagger heights of up to eight 3 liter or 5 liter bags, thus ensuring continuous flow of liquids. If you need to change an IV bag, you can do so while the other bags are still irrigating.
- Easy to clean and maintain
- Compact footprint frees up valuable space on the O.R. floor.

Affordable:

- Our IV Irrigation systems are value leaders among a few comparable models, saving you up to **\$50/procedure** on IV tubes vs. those required with far more expensive, pump driven systems.

Irrigation Systems

IVTs

Description: We offer two types of IVTs, one with 2 hooks, which holds up to 4 IV bags and one with 4 hooks, which holds up to 8 IV bags. Hooks can raise and lower from 67" to 81". They travel on black plastic twin wheel hooded swivel casters. Both IVTs come with 26 1/2" diameter inverted, 5 leg base.

Our IVT stands have the same dimensions:

Width: 12 3/4" x Height: 89" x Depth: 12 3/4"

Note: these dimensions do not include the base.



Name	IVT 2 Channels	IVT 4 Channels
Item #	IVT-2	IVT-4
Model #	0541110200	0541110100
Shipping Weight	85 lbs.*	85 lbs.*

For Custom Extender
Please email us at info@blickman.com

*Shipping weight with an asterisk is the combined shipping weight of separate parts.



Irrigation Systems

Construction Specifications

Construction Specifications

Universal Dimensions	Tower: Width: 12 3/4" x Height: 89" x Depth: 12 3/4" Base Diameter: 26 1/2" Height Range of Hooks: 67" - 81"
Materials	Hard anodized aluminum
Gauge	N/A
Welding	N/A
Finish	Stain resistant finishes
Casters	2" black plastic, twin wheel, hooded swivel casters.
Notes	

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

Note for: <i>IV Stands</i>	Parts shipped separately to avoid damage and lower costs. Some assembly required.
--------------------------------------	--

For Custom Extender

Please email us at info@blickman.com



OUR PROMISE

For over 100 years, Blickman has been known throughout the healthcare industry for exceptional product quality, and we pride ourselves on providing top notch and personalized service to ensure our customers are 100% satisfied. If you are ever not completely satisfied for any reason, we will make it right. And with our equipment backed by our Limited Lifetime Warranty, you can be certain you have made the very best choice for your healthcare facility.