

BLICKMAN



Solution Stands/ Receptacles

Roll, collect, dispose

Includes:

- Solution Stands
- Kickbucket
- Receptacle

OUR MISSION

Customer Focus

Exceed customer expectations through a commitment to continuous improvement of our products and services.

Quality Products

Offer the marketplace a broad range of high quality, cost effective, and durable products of superior craftsmanship.

Organizational Excellence

Strive for excellence in all we do by fostering a participative, team based environment, which nurtures the active involvement of all members of the organization.

Marketplace Recognition

Earn the respect, admiration, and loyalty of our people, customers, and suppliers by treating them with dignity and respect, embracing diversity, and developing and maintaining long term partnerships.

Collaborative Workplace

Encourage all employees to work safely and efficiently, contribute innovative ideas, and personally benefit from the ongoing success of the company.

The logo for Blickman, featuring the word "BLICKMAN" in a bold, white, sans-serif font. The text is contained within a dark blue rectangular box with a thin white border. The box has a slight 3D effect with a lighter blue shadow on the right side.

THE BLICKMAN DIFFERENCE

- Our high quality, heavy gauge, noise reducing, solution stands and receptacles ensure your solutions and waste products are easily collected and disposed of in these easy-to-move, easy-to-clean product.
- Our solution stands and receptacles are made out of our high quality, easy-to-clean 304 stainless steel. MR conditional options available for both. For travel and disposal, both move effortlessly and quietly on 2" swivel casters.
- Our solution stands are 33" high for compatibility with surgical table heights and are in the safe zone for reduced back stress, and our receptacles are low to the floor to pass under surgical tables and make it easy to drop waste into without reaching around the table.
- Our basins and pails are removable for both the solution stands and receptacles. They are held in place by rubber bumpers that also help prevent noise from metal on metal contact.

Key Product Features



Solution Stands

Description: We offer two types of stainless steel solution stands: a single basin model and a double basin model. Our solution stands are 33" tall offering compatibility with surgical table heights and are in the safe zone for reduced back stress. Choose between a lower shelf or an H-brace for added storage and support. MR conditional models available.



Dimensions:





Width: 15" x Height: 33" x Depth: 15"



Dimensions:

Width: 29" x Height: 33" x Depth: 15"



Name	Single Basin Solution Stand		Double Basin Solution Stand	
				
		Shelf		
Model #	7807SS	7807MR	7808SS	7808MR
Item #	0717807000	0717807066	0727808000	0727808066
Shipping Weight	20 lbs.	20 lbs.	33 lbs.	33 lbs.
	H-Brace			
Model #	7807SS-HB	7807MR-HB	7808SS-HB	7808MR-HB
Item #	0717807100	0717807166	0727808100	0727808166
Shipping Weight	19 lbs.	19 lbs.	24 lbs.	24 lbs.

Solution Basin Also Sold Separately



Name	Model #	Item #	Ship Weight
8.5 qt. Solution Basin Only	SB-1	9737370120	3 lbs.

Receptacles

Description: Description: Our stainless steel kickbucket consists of a continuous encircling rubber bumper, which helps prevent tipping, wall scuff marks, and noise, a stainless steel cross brace, which provides excellent protection as well as a secure hold of a heavy gauge 13 qt. pail, and roller blade casters for mobility.




Kickbucket
13 qt. Pail

Parts Also
Sold Separately



#KB-1



#7766NP



Receptacle
8.5 qt. Solution Basin

Parts Also
Sold Separately





#SB-1



#8766NP

KICKBUCKET

Name		Model #	Item #	Ship Weight
Lenox Kickbucket with Pail		7766SS	0817766000	12 lbs.
Lenox Kickbucket with Pail		7766MR	0817766066	12 lbs.
ORDER SEPARATELY				
Lenox Kickbucket Frame Only		7766NP	0817766001	9 lbs.
13 qt. Pail Only		KB-1	9701920000	4 lbs.

RECEPTACLE

Name		Model #	Item #	Ship Weight
Sponge Receptacle with Basin		8766SS	0828766000	11 lbs.
ORDER SEPARATELY				
Sponge Receptacle Frame Only		8766NP	0828766001	9 lbs.
8.5 qt. Solution Basin Only		SB-1	9737370120	3 lbs.



Solution Stands

Construction Specifications

Construction Specifications

Universal Dimensions	Solution Stands: 15", 29" x 33" x 15" Solution Basin: 8.5 qts. Pail: 13 qts.
Materials	304 stainless steel
Gauge	Solution Stand Frame: 1" diameter tube, 16GA Solution Stand Shelf: 20GA Pail: 22GA; Basin: 14GA
Welding	All welding by heliarc method and exposed welds are ground smooth.
Finish	Exposed surfaces are free from pit marks, weld seams, scale and are polished with #4 polish to a satin finish.
Casters	2" swivel casters.
Notes	MR conditional options available where indicated.

Shipping Information

All shipping weights are estimates and may be changed without prior notice. Specifications may be improved and dimensions changed without prior notice. Please contact customer service for estimated shipping costs at info@blickman.com.

BLICKMAN



OUR PROMISE

For over 100 years, Blickman has been known throughout the healthcare industry for exceptional product quality, and we pride ourselves on providing top notch and personalized service to ensure our customers are 100% satisfied. If you are ever not completely satisfied for any reason, we will make it right. And with our equipment backed by our Limited Lifetime Warranty, you can be certain you have made the very best choice for your healthcare facility.