

Everest & Jennings®
DURA-gel™
BASE 2G

The Dura-Gel™ BASE 2G wheelchair cushion reduces surface pressure through the use of a support layer of non-aqueous Dura-Gel™. The cushion uses a high-density base foam for durability in combination with a medium density topper.



- ❑ Features a comfortable topper with a higher density base foam
- ❑ Features a specialized Dura-Gel™ center layer surrounded by comfortable foam layers on the top and bottom of the cushion
- ❑ The cover is a 600 denier nylon with a fluid resistant backing to provide excellent durability
- ❑ Accommodates a wide range of body shapes, sizes and weights
- ❑ 250 lb maximum weight capacity
- ❑ Lightweight
- ❑ 1 year limited warranty

Latex Free

Designed to meet the requirements of HCPCS Code: E2601

8265166	16" x 16" x 2"
8265186	18" x 16" x 2"
8265188	18" x 18" x 2"
8265206	20" x 16" x 2"
8265208	20" x 18" x 2"



www.grahamfield.com

© 2012 GF Health Products, Inc.
Atlanta, Georgia 30360 GF0900045RevC12