

Everest & Jennings®

DURA-gel™

BASE 3G

The Dura-Gel™ BASE 3G wheelchair cushion reduces surface pressure through the use of a support layer of non-aqueous Dura-Gel™. The cushion uses a high-density base foam for durability in combination with a convoluted topper for comfort and airflow.



- ❑ The Dura-Gel™ BASE 3G wheelchair cushion reduces surface pressure through the use of a support layer of non-aqueous Dura-Gel™. The cushion uses a high-density base foam for durability in combination with a convoluted topper for comfort and airflow.
- ❑ Features a comfortable convoluted foam topper with a high density base foam
- ❑ Features a specialized Dura-Gel™ center layer surrounded by comfortable foam layers on the top and bottom of the cushion
- ❑ The cover has a four-way stretch top with a no-slip base fabric and has a strap at the back to secure the cushion firmly to the chair
- ❑ Accommodates a wide range of body shapes, sizes and weights
- ❑ 300 lb maximum weight capacity
- ❑ Lightweight
- ❑ 1 year limited warranty



Latex Free

Designed to meet the requirements of HCPCS Code: E2601

8565166 16" x 16" x 3"
8565186 18" x 16" x 3"
8565188 18" x 18" x 3"
8565206 20" x 16" x 3"
8565208 20" x 18" x 3"

Designed to meet the requirements of HCPCS Code: E2602

8565228 22"x18"x 3"
8565248 24"x18"x 3"



www.grahamfield.com

© 2012 GF Health Products, Inc.
 Atlanta, Georgia 30360 GF0900061RevC12