

Everest & Jennings®

DURA-gel™

SPP

The Dura-Gel™ SPP (Skin Protection & Positioning) wheelchair cushion reduces surface pressure through the use of non-aqueous Dura-Gel™ placed in three pressure zones. The cushion uses a high-density base foam in combination with a visco-elastic topper for added comfort and also positions the user with a moderate slope from front to back. Cushion also features medial thigh positioning as well as lateral positioning to provide greater stability for the user.



- ❑ Provides long-term skin protection and positioning for higher risk patients
- ❑ Reduces tissue interface pressures at critical points
- ❑ Features a specialized Dura-Gel™ center layer surrounded by comfortable foam layers on the top and bottom of the cushion
- ❑ Side wedges provide lateral stability to the lower extremities, trochanter and pelvis. Center pommel provides medial thigh support
- ❑ The soft, stretchable lightweight cover reduces shear/friction and is fluid resistant and has a strap at the back to secure the cushion firmly to the chair
- ❑ Accommodates a wide range of body shapes, sizes and weights
- ❑ 350 lb maximum weight capacity
- ❑ Lightweight
- ❑ 2 year limited warranty



Latex Free

Designed to meet the requirements of HCPCS Code: E2607

8920166	16" x 16" x 4 ½"
8920186	18" x 16" x 4 ½"
8920188	18" x 18" x 4 ½"
8920206	20" x 16" x 4 ½"
8920208	20" x 18" x 4 ½"

Designed to meet the requirements of HCPCS Code: E2608

8920228	22" x 18" x 4 ½"
8920248	24" x 18" x 4 ½"



www.grahamfield.com

© 2012 GF Health Products, Inc.
Atlanta, Georgia 30360 GF0900059RevC12