

INVACARE®
**FULL-ELECTRIC
LOW BED**



Invacare brings the safety and convenience of a low bed into the Homecare Bed Series with the **Invacare® Full-Electric Low Bed**. The split-spring design incorporated into the standard low bed profile combines the easy, one-person delivery of an Invacare bed with the benefits of a low bed, most often used when injuries from falling out of bed are a concern. The Low Bed ends are fully compatible with existing Invacare full-electric bed frame inventories and the drive shaft, simplifying inventory management for providers.

Invacare beds, mattresses and rails, excluding therapeutic support surfaces and bariatric beds, when used as a system and sold new after August 1, 2007, meet the voluntary FDA Guidelines for bed rail entrapment. Visit www.invacare.com/bedsafety for additional safety and entrapment information.

Features

High impact bed end panels are more durable, impact and scratch resistant; washable for easy cleaning

Universal low bed end can be used as either the headboard or footboard on the Low Bed

Rubber caster tread for better tracking with non-marring material

DC low voltage, compact motor system has quiet motor operation and fast cycle time on lowering height of bed for efficient bed set-up

Washable motor system and frame— Motor system with waterproof enclosures is rated IP66. Washable with standard water pressure. The bed must not be power washed

Color-coded connections on junction box for easier maintenance

Double-insulated junction box provides extra protection against electrical shock

Ergonomic hand pendant comes with durable strain relief and heavy-duty cord. Waterproof casing for easy cleaning, conveniently attaches to bed sheet with clip or bed rail with hand pendant cradle (sold separately)

Hand pendant cradle available for convenient storage of hand pendant (part no. 26975)

Interchangeability simplifies storage, reduces inventory complexity

HCPCS Code: E0297



1. Bed in lowest position
2. Bed in highest position



Invacare Corporation
www.invacare.com

USA
One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada
570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

©2012
Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and
®. All trademarks are
owned by or licensed
to Invacare Corporation
unless otherwise noted.
Specifications are subject
to change without
notification.
Form No. 04-173
Rev. 03/12
111218



Model #	Dimensions	Sleep Surface	Weight Capacity	Shipping Weight	Carton Dimensions	Regulation Compliance	Limited Warranty
5410LOW	88" L 9.5"-20" H 36" W	80" L 36" W	450 lb.			UL962 CSA Certified	5 Years Welds 2 Years Mechanical/ Electrical
5000IVC Headspring				45 lb.	36" L 44" W 8" D		
5490IVC Footspring				81 lb.	37" L 47" W 7.5" D		
5301LOW Universal Bed End				55 lb.	25" L 39" W 7" D		

All references to HCPCS codes: Providers are responsible for determining the appropriate billing codes when submitting claims to the Medicare program, and should consult an attorney or other advisor to discuss specific situations in further detail.