


PPE technical information sheet.

Product	Case Pack	Freight Specifications	Notes & Certifications
Level 1 3-Ply Disposable Surgical Masks ITEM#: MASK-3PLY	50 Pieces/Box 40 Boxes/ Case 2000 Pieces Total	Case Weight: 19lbs Case Size: 20.5" x 16.55" x 13" Cases/Skid: 24	BFE 95% FDA CE EN149 FFP2
KN95 Disposable Surgical Mask ITEM#: MASK-KN95	25 Pieces/Bag 40 Bags/ Case 1000 Pieces Total	Case Weight: 15lbs Case Size: 21" x 21" x 11" Cases/Skid :14	PFE = 99.9% @ 0.1 micron FDA CE 0121-EN149: 2001 FFP2 GB/2626-2006
Direct Splash Protection Face Shield ITEM#: FS-100PC	100 Pieces/Case	Case Size: 18.5" x 14.17" x 15.74" Case Weight: 12lbs Cases/Skid: 30	FDA Testing Report: Z87.1

3-PLY Surgical Mask



KN95 Surgical Mask



Face Shield



In addition to product packing and certification notices, this document contains certifications that have been redacted to protect proprietary sources of information.

