

# PROTEKT® 600 LIFT

## 30600-PLE | BARIATRIC POWER PATIENT LIFT

The high-performance Protekt® 600 Lift is the perfect blend of form and function for safe and easy patient transfers. This electric-powered, high quality patient transfer system provides an efficient method to safely move patients, preventing patient and caregiver injuries. Designed with the facility in mind, the lift is capable of floor lifts and fits underneath most beds in the lowest position, making it the perfect solution for transferring patients in small rooms with low beds.



(Optional Scale Item # 30300-DSC)

### SPECIFICATION FOR 30600-PLHE

weight capacity	600lb
base legs length	44.9"
base closed external width	24.9"
based open external width	42.9"
mast height	48.3"
base height	4.5"
min. lifting height	25.8"
max. lifting height	69.1"
net weight	94Lb
shipping weight	106Lb

- 600 lb. weight capacity
- Capable of lifting from the floor
- New hand-operated base opening adjusts and locks easily
- New upgraded electronic controls
- Manual emergency lowering provides safe response to power loss
- Dual front locking casters
- Clip free sling hookups improve safety and reduce downtime.
- Modular design allows easy parts replacement.
- Advanced electronics with auxiliary up/down lift controls.



### WARRANTY:

Electronics 2 years Frame 3 years

### SCALE OPTION:

Digital scale available on all Protekt® Patient Lifters. Comes with all necessary hardware.

[www.proactivemedical.com](http://www.proactivemedical.com)

