

Be Independent.™


CAR CADDIE

BECAUSE GETTING OUT OF THE CAR CAN BE EASY

PRODUCT #2080



WATCH A DEMO





Strap Adjusts Between 7"-16"



Ergonomic 5" Rubber Handle

ALSO AVAILABLE IN THE
HANDYBAR + CARCADDIE PACK
#2099



FEATURES & BENEFITS

- **Easy** – Installs and removes in seconds.
- **Universal** – Fits any framed car window.
- **Durable** – Sturdy nylon construction.
- **Versatile** – Adjustable length.
- **Secure** – No-slip rubber contour grip.

SPECIFICATIONS

Weight Capacity: 250 lbs
 Car Types: Universal Fit
 Strap Length: Adjusts from 7"–16"
 Handle Width: 5"
 Package Dimensions: 8" L x 5" W x 1" H; 0.25 lbs
 Case of 10 Dimensions: 11" L x 8" W x 6" H; 4 lbs