

EZ ADJUST BEDRAIL

BECAUSE IT'S THE WORLD'S MOST VERSATILE BED RAIL



PRODUCT #8000

















Safety Strap Holds You and the Rail in Place

Pivots 180° Down

FEATURES & BENEFITS

- Bed Rail/Hand Rail Adjustable length allows for use as both a bed rail to prevent falling out of bed or as a hand rail to asist in getting out of bed.
- Adjustable Bed rail easily adjusts from 26" to 34" to 42" after it has been placed on the bed.
- Collapsible Rail folds down to side of bed to allow user to get out
 of bed or provide space for making bed.
- Reversible Can be placed on either side of the bed
- Dual Safety Strap Secures bed rail to any home or hospital bed.
 Accommodates up to king size beds.

SPECIFICATIONS

Width of Rail: Locks in place at 26", 34" & 42" after installed on bed Height of Rail: 23" high from the base of mattress

Gap Between Rails: 3"

Bed Type: Attaches to any home or hospital bed with included safety strap

Weight of Product: 14 lbs
Weight Capacity: 400 lbs

Assembly: Installs in minutes with 4 bolts and allen wrench (included)

Package Dimensions: $28" \times 21" \times 1.5"$ Case of 4 Dimensions: $28" \times 21" \times 6"$; 60 lbs

EASY ASSEMBLY (No Additional Tools Required)



Attach L Shaped Tubes



Slide on Strap



Attach Rear Crossbar



Place on Bed