

Be Independent.™



30" SAFETY BEDRAIL

BECAUSE YOUR BED RAIL CAN PIVOT OUT OF THE WAY



PRODUCT #8050



WATCH A DEMO





Getting Out of Bed is Easier With the Hand Rail



Safety Strap Holds You and the Rail in Place



Accessorize the 30" Safety Bed Rail with this Padded Pouch #8053



Safety Bed Rail Pivots Out of Way

FEATURES & BENEFITS

- **Bed Rail/Hand Rail** - Prevents falling from bed and can be used for assistance to sit up in bed.
- **Collapsible** - Rail pivots down to side of bed to allow user to get out of bed or provide space for making bed.
- **Reversible** - Can be placed on either side of the bed.
- **Dual Safety Strap** - Secures bed rail to any sized home or hospital bed frame.

SPECIFICATIONS

Width of Rail: 30"
 Height of Rail: 23" high from the base of mattress
 Gap Between Rails: 3"
 Bed Type: Attaches to any home or hospital bed with included safety strap
 Weight of Product: 12 lbs
 Weight Capacity: 400 lbs
 Assembly: Installs in minutes with 4 bolts and allen wrench (included)
 Package Dimensions: 32" L x 21.5" W x 1.25" H 12 lbs
 Case Dimensions: 32" L x 21.5" W x 5" H; 50 lbs

EASY ASSEMBLY (No Additional Tools Required)



Attach L Shaped Tubes



Slide on Strap



Attach Rear Crossbar



Place on Bed