
Spenco RX®

DRTHOTICARCH

Orthotic Arch Supports









Firm, adjustable arch support.



Target Consumer:

Consumers with flat, low, or fallen arches. Consumers seeking firm arch support to control foot position.



Recommended For:

Treatment of arch pain. May be molded for a custom fit to support high arches.

To help reduce pronation and supination for improved stability.

Choose 3/4 length for additional toe room in tighter fitting shoes.



Best For Use In:
Casual and athletic shoes (or boots).



Components*:

Nylon, SpenCore® Material, Fabric, Arch Height: 25.80mm (1 in), Thickness at Heel: 6.29mm (1/4in)



Instructions:

Place in shoes, green side up. Trim with scissors if necessary. Adjust the shape of the orthotic by placing in boiling water for 2 minutes and then bending to desired shape.

Caution: Do not handle heated orthotic without protecting your hands. Use gloves or a towel to reshape. Allow orthotic to cool before placing in shoe.



Care:

Hand wash in warm water. Air dry.

Sizing Chart:

Size	0	1	2	3	4	5	6
U.S. Women's	3/4	5/6	7/8	9/10	11/12		
U.S. Men's			6/7	8/9	10/11	12/13	14/15
U.K.	2/3	4/5	5/6	7/8	9/10	11/12	13/14
European	32/34	36/38	38/40	40/42	42/44	44/46	46/48

Features & Benefits:

Full Length provides complete heel-to-toe comfort.

SpenCore® Material Exclusive cushioning layer absorbs shock and impact to provide maximum comfort.

Moldable Arch Insert The flexible plastic arch insert helps support and control the foot.

7 Sizes for correct arch placement for better overall fit, comfort and performance.

4-way Stretch Fabric with Silpure Antimicrobial helps prevent blisters while controlling odor.

One-Year Unconditional Guarantee

Call toll-free: 1-800-877-3626 In Canada: 1-800-387-9538



